

Effector membrane translocation biosensors reveal G protein and βarrestin coupling profiles of 100 therapeutically relevant GPCRs

Charlotte Avet^{1†}, Arturo Mancini^{2†}, Billy Breton^{2‡§}, Christian Le Gouill^{1‡}, Alexander S Hauser^{3‡}, Claire Normand², Hiroyuki Kobayashi¹, Florence Gross², Mireille Hogue¹, Viktoriya Lukasheva¹, Stéphane St-Onge¹, Marilyn Carrier¹, Madeleine Héroux¹, Sandra Morissette², Eric B Fauman⁴, Jean-Philippe Fortin⁵, Stephan Schann⁶, Xavier Leroy⁶*, David E Gloriam³*, Michel Bouvier¹*

¹Institute for Research in Immunology and Cancer (IRIC), and Department of Biochemistry and Molecular Medicine, Université de Montréal, Montréal, Canada; ²Domain Therapeutics North America, Montréal, Canada; ³Department of Drug Design and Pharmacology, University of Copenhagen, Copenhagen, Denmark; ⁴Internal Medicine Research Unit, Pfizer Worldwide Research, Development and Medical, Cambridge, United States; ⁵Pfizer Global R&D, Cambridge, United States; ⁶Domain Therapeutics, Illkirch-Strasbourg, France

Abstract The recognition that individual GPCRs can activate multiple signaling pathways has raised the possibility of developing drugs selectively targeting therapeutically relevant ones. This requires tools to determine which G proteins and β arrestins are activated by a given receptor. Here, we present a set of BRET sensors monitoring the activation of the 12 G protein subtypes based on the translocation of their effectors to the plasma membrane (EMTA). Unlike most of the existing detection systems, EMTA does not require modification of receptors or G proteins (except for G_s). EMTA was found to be suitable for the detection of constitutive activity, inverse agonism, biased signaling and polypharmacology. Profiling of 100 therapeutically relevant human GPCRs resulted in 1500 pathway-specific concentration-response curves and revealed a great diversity of coupling profiles ranging from exquisite selectivity to broad promiscuity. Overall, this work describes unique resources for studying the complexities underlying GPCR signaling and pharmacology.

Editor's evaluation

The authors describe a novel set of biosensors to assess the coupling specificity of 100 therapeutically relevant G proteins-coupled receptors (GPCRs) to various G proteins. The utility of the assay system is well-supported by the data. These tools are likely to be useful for many specific studies of individual receptors, including efforts to discover ligands that display functional selectivity (bias) between G protein pathways or between G proteins and arrestins. The work provides a rich repository of data informing on the possible effector coupling of 100 GPCRs and a set of analytical tools that could guide the development of new drugs, including efforts to discover ligands that display functional selectivity (bias) between G protein pathways or between G proteins and arrestins.

*For correspondence: xleroy@domaintherapeutics.com (XL); david.gloriam@sund.ku.dk (DEG); michel.bouvier@umontreal.ca (MB)

[†]These authors contributed equally to this work [‡]These authors also contributed equally to this work

Present address: [§]Institute for Research in Immunology and Cancer (IRIC), Université de Montréal, Montréal, Canada

Competing interest: See page 21

Funding: See page 22

Preprinted: 24 April 2020 Received: 21 September 2021 Accepted: 17 March 2022 Published: 18 March 2022

Reviewing Editor: William I Weis, Stanford University School of Medicine, United States

© Copyright Avet *et al.* This article is distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use and redistribution provided that the original author and source are

credited.

CC

Introduction

G protein-coupled receptors (GPCRs) play crucial roles in the regulation of a wide variety of physiological processes and represent one-third of clinically prescribed drugs (*Hauser et al., 2017*). Classically, GPCR-mediated signal transduction was believed to rely on linear signaling pathways whereby a given GPCR selectively activates a single G protein family, defined by the nature of its G α subunit (*Oldham and Hamm, 2008*). G α proteins are divided into four major families (G_s, G_{i/o}, G_{q/11}, and G_{12/13}) encoded by 16 human genes. Once activated, these proteins each trigger different downstream effectors yielding different biological outcomes. It has now become evident that many GPCRs can couple to more than one G protein family and that ligands can selectively promote the activation of different subsets of these pathways (*Namkung et al., 2018; Quoyer et al., 2013*). These observations extended the concept of ligand-biased signaling, which was first established for ligand-directed selectivity between β arrestin and G protein (*Azzi et al., 2003; Wei et al., 2003*), to functional selectivity between G proteins. Ligand-directed functional selectivity represents a promising avenue for GPCRs drug discovery since it offers the opportunity of activating pathways important for therapeutic efficacy while minimizing activation of pathways responsible for undesirable side effects (*Galandrin et al., 2007; Kenakin, 2019*).

To fully explore the potential of functional selectivity, it is essential to have an exhaustive description of the signaling partners that can be activated by a given receptor, providing receptor- and ligand-specific signaling signatures. Currently, few assays allow for an exhaustive pathway-specific analysis of GPCR signaling; these include BRET-based G protein activation sensors platforms (*Galés et al., 2005; Masuho et al., 2015; Maziarz et al., 2020; Mende et al., 2018; Olsen et al., 2020*) and the TGF- α shedding assay (*Inoue et al., 2019*). However, several of these platforms require modification of G protein subunits that may create functional distortions. Moreover, these assays may detect non-productive conformational rearrangements of the G protein heterotrimer as was recently reported for G₁₂ (*Okashah et al., 2020*).

Here, we describe unique sensors that do not require modification of receptors or G proteins (except for G_s) for interrogating the signaling profiles of GPCRs. The platform includes 15 pathwayselective enhanced bystander bioluminescence resonance energy transfer (ebBRET) biosensors monitoring the translocation of downstream effectors to the plasma membrane for $G_{i/o}$, $G_{q/11}$, and $G_{12/13}$, the dissociation of the G α subunit from the plasma membrane for G_s and the recruitment of βarrestin to the plasma membrane. Overall, the new ebBRET-based **E**ffector **M**embrane **T**ranslocation **A**ssays, named EMTA, provide a readily accessible large scale and comprehensive platform to study constitutive and ligand-directed GPCR signaling. The signaling signatures of 100 GPCRs using the EMTA platform also provides a rich source of information to explore the principles underlying receptor/G protein/βarrestin coupling selectivity relationships. It thus provides a unique set of tools that is complementary to previously described platforms and existing datasets, and offers a map of the coupling potentials for individual GPCR that will stimulate future studies investigating the relevance of these couplings in different physiological systems.

Results

ebBRET-based G protein effector membrane translocation assay (EMTA) allows detection of each $G\alpha$ protein subunit activation

To detect the activation of G α subtypes, we created an EMTA biosensor platform based on ebBRET (*Namkung et al., 2016; Figure 1A*). The biosensors at the heart of EMTA consist of sub-domains of the G protein-effector proteins p63-RhoGEF, Rap1GAP and PDZ-RhoGEF that selectively interact with activated G_{q/11}, G_{1/0}, or G_{12/13}, respectively. These domains were fused at their C-terminus to *Renilla* luciferase (RlucII) and co-expressed with different unmodified receptor and G α protein subtypes. Upon GPCR activation, the energy donor-fused effectors translocate to the plasma membrane to bind activated G α proteins, bringing RlucII in close proximity to the energy acceptor, *Renilla* green fluorescent protein, targeted to the plasma membrane through a CAAX motif (rGFP-CAAX), thus leading to an increase in ebBRET. The same plasma membrane translocation principle is used to measure βarrestin recruitment (*Namkung et al., 2016; Figure 1B*, top). Because no selective soluble downstream effector of G_s exists, the assay was modified taking advantage of G α_s dissociation from the plasma membrane following its activation (*Wedegaertner et al., 1996*). In this configuration,




Figure 1. EMTA ebBRET platform to monitor G protein activation and βarrestin recruitment. (**A**) Schematic of the G protein Effector Membrane Translocation Assay (GEMTA) to monitor G α protein activation. Upon receptor activation, Rlucll-tagged effector proteins (Effector-Rlucll) translocate towards and interact with active G α subunits from each G protein family, leading to increased ebBRET. (**B**) Principle of the Effector Membrane Translocation Assay (EMTA) monitoring βarrestin recruitment to the plasma membrane (*top*) and G α_s activation (*bottom*). *Top*; upon receptor activation, Rlucll-tagged βarrestins (βarrestin-Rlucll) translocate to the plasma membrane, thus increasing ebBRET with rGFP-CAAX. *Bottom*; Internalization of activated Rlucll-tagged G α_s (G α_s -Rlucll) following receptor stimulation decreases ebBRET with the membrane-anchored rGFP-CAAX.

the Rlucll is directly fused to $G\alpha_s$ (*Carr et al., 2014*). Its activation upon GPCR stimulation leads to its dissociation from the plasma membrane (*Martin and Lambert, 2016*), resulting in a reduction in ebBRET (*Figure 1B*, bottom).

The sensitivity and selectivity of the newly created G protein EMTA biosensors, were validated using prototypical GPCRs known to activate specific G α subtypes. The responses were monitored upon heterologous expression of specific G α subunits belonging to G_{i/o}, G_{q/11}, or G_{12/13} families in the absence or presence of pharmacological inhibitors and using engineered cells lacking selected G α subtypes. The dopamine D₂ receptor was used to validate the ability of the G_{i/o} binding domain of Rap1GAP (*Jordan et al., 1999*; *Meng et al., 1999*) to selectively detect G_{i/o} activation. The dopamine-promoted increase in ebBRET between Rap1GAP-RlucII and rGFP-CAAX in the presence of G $\alpha_{i/o}$ subunits was not affected by the G_{q/11}-selective inhibitor UBO-QIC (a.k.a., FR900359; *Schrage et al., 2015; Figure 2A*, left), whereas the G $\alpha_{i/o}$ family inhibitor, pertussis toxin (PTX), completely blocked the response for all members of G $\alpha_{i/o}$ family except for G α_{z} , known to be insensitive to PTX (*Casey et al., 1990; Figure 2A*, right). Gonadotropin-releasing hormone




Figure 2. Validation of EMTA ebBRET-based platform to monitor G α protein activation. (**A**) Pharmacological validation of the G $\alpha_{i/o}$ activation sensor. HEK293 cells were transfected with the D₂ receptor and the G $\alpha_{i/o}$ family-specific sensor, along with each G $\alpha_{i/o}$ subunit. Concentration-response curve using the G $\alpha_{i/o}$ activation sensor, in the presence or absence of UBO-QIC (*left*) or PTX (*right*) inhibitors. *Insets*; E_{max} values determined from

Figure 2 continued on next page

Figure 2 continued

concentration-response curves of inhibitor-pretreated cells. (B) Pharmacological validation of the $G\alpha_{\alpha'1}$ activation sensor. HEK293 cells were transfected with the GnRH receptor and the $G\alpha_{\alpha/1}$ family-specific sensor, along with each $G\alpha_{q'11}$ subunit. Concentration-response curve using $G\alpha_{q'11}$ activation sensor, in the presence or absence of UBO-QIC (*left*) or PTX (*right*) inhibitors. *Insets*; E_{max} values determined from dose-response curves of inhibitorpretreated cells. (C) Validation of the $G\alpha_{12/13}$ activation sensor. Cells were transfected with the CB₁ receptor and one of the $G\alpha_{12/13}$ activation sensors, along with the $G\alpha_{12}$ or $G\alpha_{13}$ subunits. Concentration-response curves of HEK293 cells (top) or the parental and devoid of $G_{12/13}$ ($\Delta G_{12/13}$) HEK293 cells (bottom) using the PDZ-RhoGEF-RlucII/ rGFP-CAAX (top and bottom left) or PKN-RBD-RlucII/rGFP-CAAX (bottom right) sensors, pretreated or not with UBO-QIC or PTX (top). (D) Pharmacological validation of the $G\alpha_s$ activation sensor. HEK293 cells were transfected with the GPBA receptor and the $G\alpha_s$ activation (left and central) or the EPAC (right) sensors. Left: Concentrationresponse curves using the $G\alpha_s$ activation sensor in the presence or absence of UBO-QIC or PTX, inhibitors of $G\alpha_a$ or $G\alpha_{i/o}$, respectively. Central: Concentration-response activation of the $G\alpha_s$ sensor using CTX, a $G\alpha_s$ activator. Right: Concentration-response curve using the EPAC sensor. Inset; Emax values determined from dose-response curves of inhibitors-pretreated cells. Data are expressed as BRET ratio for the concentration-response curves or expressed in % of respective control cells (E_{max} graphs) and are the mean ± SEM of 3 (A–C) or 4 (D) independent experiments performed in one replicate. Unpaired t-test (A-D): *p < 0.05 and ***p < 0.001 compared to control cells.

The online version of this article includes the following source data and figure supplement(s) for figure 2:

Source data 1. Raw data of Figure 2.

Figure supplement 1. Influence of endogenous G proteins.

Figure supplement 1—source data 1. Raw date of Figure 2—figure supplement 1.

Figure supplement 2. Validation of EMTA ebBRET-based sensors selectivity for each Ga subunit families.

Figure supplement 2—source data 1. Raw date of Figure 2—figure supplement 2.

Figure supplement 3. Influence of G protein, GPCR or effector-RlucII level expression.

Figure supplement 3—source data 1. Raw date of Figure 2—figure supplement 3.

Figure supplement 4. Kinetics of G α proteins and β arrestins recruitment promoted by the ET_A receptor.

Figure supplement 4—source data 1. Raw date of Figure 2—figure supplement 4.

Figure supplement 5. Comparison of EMTA platform and G protein activation BRET assay based on $G\alpha\beta\gamma$ dissociation.


Figure supplement 5—source data 1. Raw date of Figure 2—figure supplement 5.

Figure supplement 6. Western blots of G protein level expression in cells transfected with the EMTA ebBRET platform.

Figure supplement 6-source data 1. Original Western blot of Figure 2-figure supplement 6.


(GnRH) stimulation of the GnRH receptor (GnRHR), used as a prototypical G_{q/11}-coupled receptor, promoted ebBRET between the Rlucll-fused G_{q/11} binding domain of p63-RhoGEF (p63-RhoGEF-Rlucll; *Lutz et al., 2007; Rojas et al., 2007*) and rGFP-CAAX. The ebBRET increase observed in the presence of different G $\alpha_{q/11}$ subunits was not significantly (p = 0.077, 0.0636 and 0.073 for G_q, G₁₁, and G₁₄, respectively) affected by PTX (*Figure 2B*, right), whereas UBO-QIC completely blocked the response for all members of G $\alpha_{q/11}$ family except for G α_{15} , known to be insensitive to UBO-QIC (*Schrage et al., 2015; Figure 2B*, left). These two G protein-specific EMTA were sensitive enough to detect responses elicited by endogenous G proteins since deletion of G_{1/0} (Δ G_{1/0}) or G_{q/11} (Δ G_{q/11}) subtypes completely abolished the responses induced by D₂ or GnRHR activation in the absence of heterologously expressed G proteins (*Figure 2—figure supplement 11*). It should however be noted that relying on endogenous proteins does not allow the identification of specific members of G_{1/0} (i.e.: G₁₁, G₁₂, G₁₃, G₀₆, G₀₆, or G₂) or G_{q/11} (i.e.: G_q, G₁₁, G₁₄, or G₁₅) families.

The selectivity of the G_{12/13} binding domain of PDZ-RhoGEF (*Fukuhara et al., 2001*) was confirmed using the cannabinoid receptor type 1 (CB₁). The ebBRET between PDZ-RhoGEF-RlucII and rGFP-CAAX in the presence of G α_{12} or G α_{13} promoted by the cannabinoid agonist WIN-55,212–2 was not affected by UBO-QIC (*Figure 2C*, top left), nor PTX (*Figure 2C*, top right). Given the lack of selective G_{12/13} pharmacological inhibitor, we used HEK293 cells genetically deleted for G α_{12} and G α_{13} proteins (Δ G_{12/13}) to further confirm the response selectivity. As expected, PDZ-RhoGEF-RlucII/rGFP-CAAX ebBRET was observed only following reintroduction of either G α_{12} (Δ G_{12/13}+G₁₂) or G α_{13} (Δ G_{12/13}+G₁₃) (*Figure 2C*, bottom left).



Video 1. BRET-based imagery of p63-RhoGEF-RlucII recruitment to the plasma membrane upon AT₁ activation. HEK293 cells expressing the p63-RhoGEF-RlucII/rGFP-CAAX sensors with $G\alpha_q$ and AT₁ were stimulated with Angiotensin II. BRET levels (the ratio of the acceptor photon count to the total photon count) are expressed as a color code (lowest being black and purple, and highest being red and white). https://elifesciences.org/articles/74101/figures#video1

Biochemistry and Chemical Biology




Video 2. BRET-based imagery of Rap1GAP-RlucII recruitment to the plasma membrane upon D_2 activation. HEK293 cells expressing the Rap1GAP-RlucII/rGFP-CAAX sensors with $G\alpha_{i2}$ and D_2 were stimulated with dopamine. BRET levels (the ratio of the acceptor photon count to the total photon count) are expressed as a color code (lowest being black and purple, and highest being red and white). https://elifesciences.org/articles/74101/figures#video2

The $G_{12/13}$ coupling of CB₁ was further confirmed by monitoring the recruitment of PKN to the plasma membrane (*Figure 2C*, bottom right) in agreement with previous reports (*Inoue et al., 2019*).

To further assess the selectivity of each EMTA biosensor, we took advantage of the fact that the

endothelin-1 receptor (ET_A) can activate $G_{q/11}$, $G_{i/o}$, and G_{12/13} family members. As shown in Figure 2figure supplement 2, only over-expression of the Ga family members corresponding to their selective effectors (Rap1GAP for G_{i/o}, p63-RhoGEF for $G_{q/11}$, and PDZ-RhoGEF for $G_{12/13}$) significantly increased the recruitment of the effector-RlucII to the plasma membrane. A recent study (Chandan et al., 2021) showed that G_{i/o} can also activate full length PDZ-RhoGEF. Although the domain of PDZ-RhoGEF required for this activation has not been identified yet, the selectivity of our PDZ-RhoGEF sensor for $G_{12/13}$ vs. all other G protein families most likely results from the fact that we used a truncated version of PDZ-RhoGEF that only contains the $G_{12/13}$ binding domain and lacks the PDZ domain involved in protein-protein interaction, the actin-binding domain and the DH/PH domains involved in GEF activity and RhoA activation (Aittaleb et al., 2010).

It should be noted that in the heterologous expression configuration, competition with endogenous G proteins did not occur to a significant extent since the potencies of the responses to a given G protein subtype were not affected by genetic deletion of the different G protein



Video 3. BRET-based imagery of PDZ-RhoGEF-RlucII recruitment to the plasma membrane upon TP α R activation. HEK293 cells expressing the PDZ-RhoGEF-RlucII/rGFP-CAAX + G α_{13} and TP α R were stimulated with U46619. BRET levels (the ratio of the acceptor photon count to the total photon count) are expressed as a color code (lowest being black and purple, and highest being red and white).

https://elifesciences.org/articles/74101/figures#video3

family members (Figure 2—figure supplement 1 and Supplementary file 1A). Similarly, overexpression of G proteins, GPCRs or effectors-RlucII did not affect the potencies of the responses observed (Figure 2—figure supplement 3 and Supplementary file 1B-D), indicating that, in our experimental conditions, overexpression of the different components of EMTA sensors must likely not bias the coupling response. In addition to spectrometric assessment of coupling selectivity (above) and activation kinetics (Figure 2—figure supplement 4), EMTA allows to image the real-time recruitment of the G protein effectors to the plasma membrane (Videos 1–3) thus providing spatiotemporal resolution for the imaging detection of $G\alpha_{i/o}$, $G\alpha_{a/11}$, and $G\alpha_{12/13}$ activation.

The sensitivity of the EMTA platform is illustrated by a direct side-by-side comparison of the signals detected with EMTA vs. BRET assays based on $G\alpha\beta\gamma$ dissociation ($G\alpha\beta\gamma$) (*Galés et al., 2005*; *Galés et al., 2006*; *Olsen et al., 2020*), that reveals a significantly larger assay windows for EMTA for the 6 $G\alpha$ subunits tested for eight selected receptors, (*Figure 2—figure supplement 5*).

For the G α_s translocation biosensor, the bile acid receptor (GPBA) was chosen for validation (*Kawa-mata et al., 2003*). As expected, lithocholic acid stimulation resulted in a concentration-dependent decrease in ebBRET between G α_s -Rlucll and rGFP-CAAX (*Figure 2D*, left). Cholera toxin (CTX), which directly activates G α_s (*De Haan and Hirst, 2004*), led to a decrease in ebBRET (*Figure 2D*, center), confirming that loss of G α_s plasma membrane localization results from its activation. The potency of lithocholic acid to promote G_s dissociation from the plasma membrane was well in line with its potency to increase cAMP production as assessed using a BRET²-based EPAC biosensor (*Leduc et al., 2009*; *Figure 2D*, right). The G_s-plasma membrane dissociation ebBRET signal was not affected by UBO-QIC or PTX (*Figure 2D*, left), confirming the selectivity of the biosensor.

Signaling signatures of one hundred therapeutically relevant receptors reveals distinct G protein and β arrestin selectivity profiles

We used EMTA to assess the signaling signature of a panel of 100 human GPCRs that are either already the target of clinically used drugs (74 receptors), considered for pre- or clinical drug development (6 receptors), or pathophysiologically relevant (*Supplementary file 2A*). To establish the coupling potentials for each receptor, we quantified its ability to activate 15 pathways: $G\alpha_s$, $G\alpha_{i1}$, $G\alpha_{i2}$, $G\alpha_{oA}$, $G\alpha_{oB}$, $G\alpha_2$, $G\alpha_{12}$, $G\alpha_{13}$, $G\alpha_q$, $G\alpha_{11}$, $G\alpha_{15}$ and β arrestin 2 as well as β arrestin 1 and 2 in the presence of GRK2 (*Supplementary file 3*). E_{max} and pEC₅₀ values were determined (*Supplementary file 2*) and, based on the pre-determined threshold criteria (Emax ≥mean of vehicle-stimulated +2*SD; see Materials and methods), a 'yes or no' agonist-dependent activation was assigned to each signaling pathway and summarized using radial graph representations (*Figure 3—figure supplement 1*). To assess whether endogenous receptors could contribute to the observed responses, assays were also carried out in cells not transfected with the studied receptor (*Figure 3—figure supplement 2*). When an agonist-promoted response was observed in non-transfected parental HEK293 cells, this response was not considered as a receptor-specific response (see Materials and methods).

To compare the signaling profiles across all receptors and pathways and to overcome differences in receptor expression levels and individual biosensor dynamic windows, we first min-max normalized E_{max} and pEC₅₀ values (between 0 and 1) across receptors as a function of a reference receptor yielding the largest response for a given pathway (Figure 3A, left). Then, these values were again min-max normalized (between 0 and 1) for the same receptor across pathways, using the pathway with the largest response for this receptor as the reference (Figure 3A, right; see description in Materials and methods). Such double normalization allows direct comparison of the coupling efficiency to different G proteins for a given receptor and across receptors for a given G protein. This coupling efficiency is summarized as heatmaps (Figure 3B) that reveals a high diversity of signaling profiles. The selectivity toward the different G protein families varies considerably among GPCRs (Figure 4). In our dataset, which is the first using unmodified GPCRs and G α proteins (except for G_s), 29% of the receptors coupled to only one family, whereas others displayed more promiscuity by coupling to 2, 3, or 4 families (36%, 25%, and 10%, respectively). Receptors coupling to a single G protein family favored the members of the $G_{i/o}$ family. Indeed, 27% of the receptors coupling to $G_{i/o}$ only activated this subtype family in comparison to 0, 2.4 and 9.1% for receptors activating $G_{12/13}$, $G_{\alpha/11}$, and G_s , respectively, thus displaying more promiscuous coupling. A detailed comparative analysis of the selectivity profiles that we observed using the EMTA sensors with that of the chimeric G protein-based assay developed by Inoue et al., 2019 and the IUPHAR/BPS Guide to Pharmacology database (GtP;




Figure 3. Heatmaps illustrating the diversity of receptor-specific signaling signatures detected with the EMTA ebBRET platform. (**A**) First, values within each pathway were normalized relative to the maximal response observed across all receptors (max = 1; *left*). These values were then normalized across pathways for the same receptor, with the highest-ranking pathway serving as the reference (max = 1; *right*). (**B**) Heatmap representation of double

Figure 3 continued on next page

Figure 3 continued

normalized E_{max} (*left*) and pEC₅₀ (*right*) data. Empty cells (grey) indicate no detected coupling. IUPHAR receptor names are displayed.

The online version of this article includes the following source data and figure supplement(s) for figure 3:

Figure supplement 1. Receptor-specific signaling signatures.

Figure supplement 2. Detection of endogenous receptor-mediated responses with the EMTA ebBRET platform in HEK293 cells.

Figure supplement 2—source data 1. Raw data of Figure 3—figure supplement 2.

Figure supplement 3. Validation of G_{12/13} and G₁₅ signaling for the newly characterized GPCRs.

Figure supplement 3—source data 1. Raw data of Figure 3—figure supplement 3.

https://www.guidetopharmacology.org/) is presented in the accompanying paper (*Hauser et al., 2022*). *Supplementary file 2C* allows a direct comparison of the relative potency determined using EMTA for both the new and the already known (i.e.: identified in GtP database) couplings. As can be seen in the table, although in many cases the potency for the novel couplings is lower, this is not a universal finding since for some receptors, the pEC₅₀s for the new couplings are similar (ex: G_{12} for CB₁; G_{13} for serotonin 5-HT_{2C}; $G_{12/13}$ for adenosine 2A (A_{2A}) and prostaglandin E1 (EP₁) receptors; $G_{i/o}$ for corticotropin-releasing hormone receptor 1 (CRFR1), ET_A and G protein-coupled receptor 39 (GPR39))




Figure 4. The EMTA ebBRET platform has a unique ability to uncover coupling selectivity between G protein families. (A) Venn diagram showing the numbers of receptors coupled to each G protein family in the EMTA ebBRET biosensor assay. (B) Evaluation of receptors coupling promiscuity: number of receptors that couple to members of 1, 2, 3, or 4 G protein families. (C) Determination of G protein subunit coupling frequency: number of receptors that activate each G α subunit. (D) Proportion of receptors recruiting βarrestins: number of receptors that do not recruit (-/-) or that recruit either (+/- or -/+) or both (+/+) βarrestin isotypes. All data are based on double normalized E_{max} values from *Figure 3*.

The online version of this article includes the following source data and figure supplement(s) for figure 4:

Source data 1. Raw data of Figure 4.

Figure supplement 1. G protein subtypes distribution across the 100 GPCRs profiled with the EMTA ebBRETbased platform.

Figure supplement 1—source data 1. Raw data of Figure 4—figure supplement 1.

or higher (ex: G_z for serotonin 5-HT₂₈; G_{15} for adenosine 3 (A₃) and melanocortin 3 (MC3R) receptors; G_{12} for bradykinin 2 (B₂), cholecystokinin A (CCK₁), chemokine receptor 6 (CCR6) and ET_A receptors; $G_{12/13}$ for CRFR1 and GPR68) than those for the canonical ones. Interestingly, in many instances the potency for the newly uncovered couplings are similar to those for βarrestins, which is generally lower than for their canonical G proteins, a finding consistent with the role of βarrestins in signaling arrest at the plasma membrane. The potency differences observed for the activation of different G protein subtypes by a given receptor may lead to preferential activation of some pathways over others. This relative selectivity is likely to be influenced by tissue-dependent G protein subtype expression levels. The physiological consequences of such selectivity remain to be investigated.

When examining the frequency of coupling for each G α subunit family (Figure 4C), the G_{i/o} family members were the most commonly activated, with 89% of the tested receptors activating a Gi/o family member. In contrast, only 33%, 49%, and 45% of the receptors activate G_s , $G_{12/13}$, or $G_{q/11}$ (excluding $G\alpha_{15}$ family members, respectively. Not surprisingly, and consistent with its reported promiscuous coupling, $G\alpha_{15}$ was found to be activated by 81% of the receptors. For some receptors, we also observed preferential coupling of distinct members within a subtype family (Figure 3-figure supplement 1). For instance, 33% of Gi/o-coupled receptors can couple to only a subpopulation of the family (Figure 4—figure supplement 1A). For the $G_{a/11}$ family, only 44% activate all family members with 45% activating only $G\alpha_{15}$ and 11% engaging only two or three members of the family. A matrix expressing the % of receptors engaging a specific G α subtype that also activated another subtype, is illustrated in Figure 4—figure supplement 1B. When considering individual families, considerable variation within the $G_{i/o}$ family was observed. The greatest similarities were observed between $G\alpha_{oB}$ and either $G\alpha_{oA}$ or $G\alpha_z$, and the lowest between $G\alpha_1$ and $G\alpha_z$. A striking example of intra-family coupling selectivity is the serotonin 5-HT_{2B} that activates only $G\alpha_{oB}$ and $G\alpha_{z}$ and GPR65 that selectively activates $G\alpha_{oB}$. Similarly, when considering the ligand-promoted responses above our threshold criteria (see Materials and methods), histamine H_2 and MC3R receptors show preferred activation of $G\alpha_{_{OB}}$ and $G\alpha_{_{zr}}$, whereas the prostaglandin F (FP) and neuropeptide Y5 (Y₅) receptors preferentially activate $G\alpha_{oB}$, $G\alpha_{oA}$, and $G\alpha_{z}$. Even when all members of a given family are found to be activated, some receptors activate specific family members with greater potencies (Supplementary file 2C).

When considering β arrestin recruitment, our analysis shows that 22% of receptors did not recruit β arrestin 1 or 2, even in the presence of overexpressed GRK2 (*Figure 4D*). Among the receptors able to recruit β arrestins, only a very small number selectively recruited β arrestin1 (1.3%) or β arrestin2 (6.4%), most of them recruiting both β arrestins in the presence of GRK2 (92.3%) (*Figure 4D*). Overexpression of GRK2 potentiated the recruitment of β arrestin2 for 68% of receptors highlighting the importance of GRK2 expression level in determining β arrestin activation (*Supplementary files 3 and 2*).

Comparison with previous datasets reveals commonalities and crucial differences

We compared the signaling profiles obtained here with those presented by **Inoue et al.**, **2019** and the GtP dataset. Of note, this comparison only considers the final reported couplings that in the Inoue's study were based on the criteria of positive coupling if LogRAi ≥ -1 and negative coupling if LogRAi ≤ -1 , and is influenced by the different cut-offs and normalization used in the two studies. A comparison of couplings using common E_{max} standard deviation cut-off, quantitative normalization and aggregation of G proteins into families is provided in the accompanying paper (**Hauser et al., 2022**). As can be seen in **Supplementary file 4A**, among the 70 receptors common to both studies, less couplings were detected in our study than reported in Inoue et al. for G α_s (21 vs. 28), G α_1 (54 vs. 56), G α_q (31 vs. 34), and G α_{14} (36 vs. 40). In contrast, more receptors activating G α_{12} (29 vs. 23), G α_o (59 vs. 41), G α_{13} (30 vs. 15), G α_z (52 vs. 37), and G α_{15} (62 vs. 15) were detected in our study. When comparing with data collected in GtP, that reports couplings grouped for G protein families (*i.e.*: G_s, G_{i/o}, G_{q/11}, or G_{12/13}) and not at the single G protein subtype level, we detected less couplings than what was reported in GtP for G α_s (32 vs. 37), but more for G $\alpha_{i/o}$ (89 vs. 69), G $\alpha_{q/11}$ (81 vs. 48), and G $\alpha_{12/13}$ (47 vs. 10), among the 99 receptors common to both datasets (**Supplementary file 4B**).

Altogether, the comparative analysis reveals 64% and 69% identity of couplings between the EMTA and Inoue's or GtP datasets, respectively. Each dataset reporting unique couplings and missing couplings found in the other two datasets. The reasons for these differences are plausibly due to

intrinsic differences in the assays used. For instance, for $G_{12/13}$ and G_{15} specifically, the difference with the GtP dataset most likely results from the fact that in most cases $G_{12/13}$ or G_{15} activation were determined indirectly since, until their recent description ($G_{12/13}$: **Quoyer et al., 2013**; **Schrage et al., 2015**; G_{15} :**Inoue et al., 2019**; **Olsen et al., 2020**), no robust readily available assay existed to monitor the activation of these G proteins.

Validation of newly identified G_{12/13} and G₁₅ couplings

Given the overrepresentation of both $G_{12/13}$ and G_{15} couplings, obtained with the EMTA assays vs. those reported by Inoue et al. and the GtP datasets, the validity of the EMTA assay to detect real productive couplings, was confirmed using orthogonal assays for selective examples not reported in the two other datasets. For $G_{12/13}$, we used the PKN-based BRET biosensor detecting Rho activation downstream of either $G_{12/13}$ or $G_{q/11}$ (*Namkung et al., 2018*) and the MyrPB-Ezrin-based BRET biosensor detecting the activation of Ezrin downstream of $G_{12/13}$ (*Leguay et al., 2021*), both in the absence of heterologously expressed G proteins. Ligand stimulation of FP and CysLT₂ receptors led to Rho and ezrin activation (*Figure 3—figure supplement 3A*), that were insensitive to the $G_{q/11}$ inhibitor YM-254890, confirming that these receptors activate $G\alpha_{12/13}$.

For newly identified G_{15} couplings, we took advantage of the lack of $G\alpha_{15}$ in HEK293 cells and assessed the impact of $G\alpha_{15}$ heterologous expression on receptor-mediated calcium responses (*Figure 3—figure supplement 3B*). For prostaglandin E2 (EP₂) and κ -opioid (κ OR) receptors, which couple to G_{15} but no other $G_{q/11}$ members, expression of $G\alpha_{15}$ significantly increased the PGE2- and Dynorphin A- promoted calcium responses. For α_{2A} adrenergic ($\alpha_{2A}AR$) and vasopressin 2 (V₂) receptors that couple other $G_{q/11}$ family members, treatment with YM-254890 completely abolished the agonistpromoted calcium response in the absence of $G\alpha_{15}$. In contrast, the calcium response evoked by $\alpha_{2A}AR$ and V₂ agonists following $G\alpha_{15}$ expression was completely insensitive to YM-254890 (*Figure 3—figure supplement 3B*), confirming that these receptors can activate this YM-254890-insensitive G protein subtype (*Takasaki et al., 2004*).

EMTA platform detects constitutive receptor activity and biased signaling

We went on to assess the ability of the EMTA platform to detect receptor constitutive activity. Transfection of increasing amounts of adenosine A_1 receptor (A_1) led to a receptor-dependent increase in basal ebBRET of the $G\alpha_{i2}$ -activation sensor (*Figure 5A*, left), reflecting A_1 constitutive activity. The A_1 inverse agonist DPCPX (*Lu et al., 2014*) dose-dependently decreased the constitutive A_1 -mediated activation of $G\alpha_{i2}$ (*Figure 5A*, left), indicating that EMTA can detect inverse agonism. Although we can not exclude that the high basal activity was observed for A_1 but not A_3 , despite a similar potency of adenosine to activate these two receptors subtypes (see *Figure 5—figure supplement 1A*), supports the notion that the increased basal activity reflects A_1 constitutive activity.

To further confirm that the platform can adequately detect inverse agonism, a second receptor for which no endogenous ligand should be present in the media, the CB₁ receptor, was used. As illustrated in *Figure 5A* (right), increase CB₁ expression led to a ligand-independent constitutive activation of G_z, that could be completely blocked by the CB₁ inverse agonist rimonabant.

EMTA also faithfully detected biased signaling. Indeed, as previously reported (**Namkung et al., 2018**; **Wei et al., 2003**), angiotensin analogs such as SII, saralasin or TRV027 displayed biased signaling by promoting efficient β arrestin2 recruitment but marginal or no $G\alpha_q$, $G\alpha_{i2}$, or $G\alpha_{13}$ activation as compared to angiotensin II that activated all G proteins and β arrestin2 (*Figure 5B*). The platform was also used to identify biased signaling resulting from single nucleotide polymorphisms. As shown in *Figure 5C*, two naturally occurring variants of human GPR17 (isoform 2) localised in the TM3 E/ DRY motif resulted in altered functional selectivity profiles. Whereas the Asp128Asn variant displayed WT-like activity on $G\alpha_{i2}$, it lost the ability to activate $G\alpha_q$ and β arrestin2. In contrast, variant Arg129His at the neighboring position resulted in an increased constitutive β arrestin2 recruitment and a loss of $G\alpha_{i2}$ and $G\alpha_q$ protein signaling.




Figure 5. Multiple applications using the EMTA ebBRET platform. (A) Inverse agonist activity detection. Left: Ga₂ activation in HEK293 cells transfected with the Rap1GAP-RlucII/rGFP-CAAX sensors with untagged G α_{12} and increasing amount of A₁ receptor plasmid. Data are expressed in % of response obtained in control cells (0 ng of A1) and are the mean ± SEM of 4-6 independent experiments performed in two replicates. One Way ANOVA test: ***p < 0.001 compared to control cells. HEK293 cells expressing the G α_2 activation sensor and control (Mock) or A₁ receptor plasmid were stimulated (10 min) with increasing concentrations of the indicated compound. Data are expressed in % of constitutive response obtained in vehicle-treated A1 transfected cells and are the mean ± SEM of 4-6 independent experiments performed in one replicate. Right: Gaz activation in HEK293 cells transfected with the Rap1GAP-RlucII/rGFP-CAAX sensors with untagged G α_z and increasing amount of CB₁ receptor plasmid. Data are expressed in % of response obtained in control cells (0 ng of CB1) and are the mean ± SEM of 4 independent experiments performed in one replicate. One Way ANOVA test: ***p < 0.001 compared to control cells. HEK293 cells expressing the G α_z activation sensor and increasing amount of CB₁ receptor plasmid were directly stimulated (10 min) with increasing concentrations of the CB1 inverse agonist rimonabant. Data are expressed as % of the response obtained in control cells (0 ng of CB₁) treated with vehicle and are the mean ± SEM of 4 independent experiments performed in one replicate. (B) Ligand-biased detection. Concentration-response curves of AT₁ for the endogenous ligand (Angiotensin II, AngII) and biased agonists [Sar1-Ile4-Ile8] AngII (SII), saralasin or TRV027. G protein and βarrestin2 signaling activity were assessed by EMTA platform. Data are expressed in % of maximal response elicited by Angll and are the mean ± SEM of 3-6 independent experiments performed in one replicate. (C) Functional selectivity of naturally occurring receptor variants. Concentration-response curves for WT or E/DRY motif Asp128Asn and Arg129His variants of GPR17 upon agonist stimulation in HEK293 cells co-expressing the indicated EMTA biosensor. Data are expressed in % of maximal response elicited by WT receptor and are the mean ± SEM of 3 independent experiments performed in one replicate.

The online version of this article includes the following source data and figure supplement(s) for figure 5:

Source data 1. Raw data of Figure 5.

Figure supplement 1. Modulation of ligand-promoted response detected by EMTA ebBRET platform by receptor constitutive activity.

Figure supplement 1—source data 1. Raw data of Figure 5—figure supplement 1.

Combining $G_{\scriptscriptstyle Z}$ and $G_{\scriptscriptstyle 15}$ biosensors for safety panels and systems pharmacology

The G protein coupling profiles obtained for the 100 GPCRs revealed that 95% of receptors activate either $G\alpha_z$ (73%) or $G\alpha_{15}$ (81%). Measuring activation of both pathways simultaneously provides an almost universal sensor applicable to screening. Combining the two sensors (Rap1GAP-Rlucll/ p63-RhoGEF-Rlucll/rGFP-CAAX) in the same cells allowed to detect ligand concentration-dependent activation of a safety panel of 24 GPCRs, that are well established as contributors to clinical adverse




Figure 6. Detection of direct and indirect (*trans*) mechanisms of ligand polypharmacology using the G_z/G_{15} biosensor. (**A**) Test of the G_z/G_{15} biosensor on a safety target panel. ebBRET signal was measured before and after stimulation with the indicated ligand in HEK293 cells transfected with the combined G_z/G_{15} biosensor and one of the 24 receptors listed. (**B**) Cross-activation of D_2 and $\alpha_{2A}AR$ by other natural ligands. For the agonist mode read, HEK293 cells expressing D_2 or $\alpha_{2A}AR$ and either the $G\alpha_{i2}$, $G\alpha_{oB}$, or the βarrestin2 + GRK2 sensors were stimulated with increasing concentrations of the indicated ligand. For the antagonist mode read, cells were pretreated with increasing concentrations of the selective D_2 antagonist eticlopride or the selective $\alpha_{2A}AR$ antagonist WB4101 before stimulation with an EC_{80} of the indicated ligand. Data are the mean ± SEM from 3-4 independent experiments performed in one replicate and expressed in % of the response elicited by dopamine or noradrenaline for D_2 and $\alpha_{2A}AR$ expressing cells, respectively. (**C**) Indirect (*trans*) activation of CB₁ by acetylcholine. For the agonist mode read, HEK293 cells expressing CB₁ and the Rap1GAP-RlucII/ rGFP-CAAX sensors with untagged $G\alpha_{oB}$ were stimulated with increasing concentrations of the indicated ligand. For the antagonist mode read, same cells were pretreated or not with increasing concentrations of the CB inverse agonist AM-630 (*left*) or the cholinergic antagonist atropine (*central*) before stimulation with an EC₈₀ of the indicated ligand (*right*). Data are the mean ± SEM from 3-5 independent experiments performed in one replicate and expressed in % of the response elicited by WIN55,212–2.

The online version of this article includes the following source data and figure supplement(s) for figure 6:

Source data 1. Raw data of Figure 6.

Figure supplement 1. Combined G_z/G₁₅ biosensor.

Figure supplement 1—source data 1. Raw data of Figure 6—figure supplement 1.

Figure supplement 2. Validation of direct activation of $\alpha_{2A}AR$ by dopamine.

Figure supplement 2—source data 1. Raw data of Figure 6—figure supplement 2.

drug reactions (Bowes et al., 2012; Figure 6-figure supplement 1). Indeed, the G₂/G₁₅ sensor captured the activation of receptors largely or uniquely coupled to either $G\alpha_z$ (e.g. CB_2) or $G\alpha_{15}$ (e.g. A_{2A} and A_{2B}), as well as receptors coupled (to varying degrees) to both pathways. The usefulness of the G_z/G_{15} combined sensor to detect off-target ligand activity is illustrated in *Figure 6A*. Most ligands tested were specific for their primary target(s). However, certain ligands displayed functional cross-reactivity with GPCRs other than their cognate targets. These included the activation of the $\alpha_{2A}AR$ by dopamine and serotonin, the D₂ by noradrenaline and serotonin, and of the CB₁ and CB₂ receptors by acetylcholine (Figure 6B-C). The activation of D₂ by noradrenaline and serotonin was confirmed by the ability of the D2-family selective antagonist eticlopride to block the dopamine-, serotonin-, and noradrenaline-promoted responses detected using the combined G_{ν}/G_{15} or the G_{ν} and G_{oB} -selective sensors and β arrestin2 sensor (*Figure 6B*, top). Similarly, use of the α_2AR selective antagonist, WB4101, allowed to confirm that dopamine can activate $G\alpha_{_{12}}$, $G\alpha_{_{OB}}$ and β arrestin2 through the $\alpha_{2A}AR$ (Figure 6B, bottom). Such pleiotropic activation of different monoaminergic receptors by catecholamines and serotonin has been previously observed (Roth et al., 2004; Sánchez-Soto et al., **2016**; Sunahara et al., 1991). Direct activation of the $\alpha_{2A}AR$ by dopamine was confirmed by showing that treatment with the D₂-family receptor selective antagonist eticlopride had negligible effect on dopamine-mediated activation of $G\alpha_{12}$ and $G\alpha_{0B}$ in cells heterologously expressing $\alpha_{2A}AR$, confirming that the response did not result from the activation of endogenously expressed dopamine receptor. In contrast, eticlopride blocked the activation of $G\alpha_{i2}$ and $G\alpha_{oB}$ in cells heterologously expressing D_2 (Figure 6—figure supplement 2).

These cross-reactivity may be direct (i.e. via direct binding of a ligand to its non-cognate receptor) as suggested above, or indirect (e.g. 'trans', via ligand activation of its canonical receptor, leading to subsequent secretion of factors that activate the non-canonical target). One such example of trans-activation is provided by the activation of cannabinoid CB1 and CB2 receptors by acetylcholine (detected by the $G_{z/15}$ and confirmed with the G_{oB} sensors; Figure 6A and C). Indeed, the activation was completely inhibited by both the CB inverse agonist AM-630 and by the cholinergic antagonist atropine (Figure 6C, left). Yet the response evoked by the CB selective agonist WIN55,212 2 was not blocked by atropine (Figure 6C, center). $G\alpha_{oB}$ activation by acetylcholine did not result from direct activation of endogenous muscarinic receptors since no $G\alpha_{oB}$ response was observed in parental cells (Figure 3-figure supplement 2). Given that the M₃ muscarinic receptor, which is endogenously expressed at relatively high levels in HEK293 cells (Atwood et al., 2011), is strongly coupled to the $G_{a/11}$, CB₁-expressing cells were pretreated with $G_{a/11/14}$ inhibitor UBO-QIC prior to stimulation with acetylcholine. UBO-QIC pre-treatment blocked acetylcholine- but not WIN55,212–2-mediated $G\alpha_{OB}$ activation (Figure 6C, right). These results demonstrate that CB1 activation by acetylcholine is indirect and potentially involves the secretion of an endogenous CBR ligand following activation of $G_{\alpha/11}$ by endogenous muscarinic acetylcholine receptors. The combined Gz/G15 sensor is therefore a useful tool to identify interplay between receptors and to explore systems pharmacology resulting from such cross-talks.

Discussion

This study describes the development and validation of a genetically encoded ebBRET-based biosensor platform allowing live-cell mapping of GPCR-G protein coupling preferences covering 12 heterotrimeric G proteins. The novel EMTA biosensors were combined with previously described ebBRET-based β arrestin trafficking sensors (*Namkung et al., 2016*), providing an unprecedented description of GPCR signaling partner couplings. In addition to providing a resource to study GPCR functional selectivity (*Pándy-Szekeres et al., 2022*), the sensors provide versatile and readily usable tools to study, on a large-scale, pharmacological processes such as constitutive activity, inverse agonism, ligand-biased signaling, and signaling cross-talk.

Our EMTA-based biosensor platform offers several advantages relative to other available approaches. First, EMTA provides direct real-time measurement of proximal signaling events following GPCR activation (i.e. $G\alpha$ protein activation and β arrestin recruitment) and resulting in lower level of amplification than those of assays relying on enzymatic activity of downstream effectors (i.e.: adenylyl cyclase or phospholipase C) or artificial detection systems (i.e.: gene-reporter or TGF- α shedding assays) that measure signal accumulation sometimes following extended incubation times. In addition,




Figure 7. Detection of endogenous receptor- and/or G protein-mediated responses in cells with the EMTA ebBRET platform. Concentration-dependent activation of $G\alpha_{i2}$ protein by (**A**) endogenous S1P₁ receptor in iPSC-derived cardiomyocytes transfected with heterologous $G\alpha_{i2}$, (**B**) endogenous FPR2 in promyelocytic HL-60 cells transfected with heterologous $G\alpha_{i2}$, (**C**) endogenous FPR2 in promyelocytic HL-60 cells with endogenous $G_{i/o}$ proteins and (**D**) endogenous PAR2 receptor in HEK293 cells with endogenous $G_{i/o}$ proteins. In all cases, cells were co-transfected with the Rap1GAP-Rlucll/rGFP-CAAX biosensor. Data are the mean ± SEM of 3-4 independent experiments performed in one replicate and are expressed as BRET² ratio in percentage of response induced by vehicle.

The online version of this article includes the following source data for figure 7:

Source data 1. Raw data of Figure 7.

measuring proximal activity reduces the risk of cross-talks between pathways that may complicate data interpretation when considering downstream signaling as the readout (*Mancini et al., 2015*).

Second, EMTA uses native untagged GPCRs and G protein subunits (except for G_s), contrary to protein complementation (Laschet et al., 2019), FRET/BRET-based Gαβy dissociation/receptor-G protein interaction (Bünemann et al., 2003; Galés et al., 2005; Galés et al., 2006; Hoffmann et al., 2005; Namkung et al., 2018; Olsen et al., 2020) or TGF-α shedding (Inoue et al., 2019) assays. Modifying these core-signaling components could alter responses, complicate interpretation and explain some of the discrepancies observed between the EMTA platform and other approaches used to study G protein activation. Moreover, the ability to work with unmodified receptors and G proteins (except for G_e) offers numerous advantages. First, it allows for the detection of endogenous GPCR signaling in either generic HEK293 cells (Figure 3-figure supplement 2) or more physiologically relevant cell lines such as induced pluripotent stem cell (iPSC)-derived cardiomyocytes (Figure 7A) and promyelocytic HL-60 cells (Figure 7B). Further it allows, in cells expressing sufficient endogenous level of the G proteins of interest, to detect activation of both native receptor and G proteins with no need of overexpression (Figure 7C-D). This is illustrated by the ability to detect the recruitment of Rap1GAP upon activation of the endogenous $G_{i/a}$ family members by the formyl peptide receptor 2 (FPR2) in HL-60 cells (Figure 7C) or protease-activated receptor-2 (PAR2) in HEK293 cells (Figure 7D). The ability to detect the activation of endogenous G protein was also illustrated in Figure 2-figure supplement 11, where the responses elicited by agonist stimulation were lost in cells genetically deleted of the G protein engaged by the studied receptor (i.e.: $G_{\alpha'11}$ or $G_{i/\alpha}$ families). Recently, another BRET-based approach (Maziarz et al., 2020), taking advantage of a synthetic peptide recognizing the GTP-bound form of G α subunits, also allows the detection of native G protein activation, offering alternative means to probe coupling selectivity profiles for both endogenously and heterologously expressed GPCRs.

Finally, similarly to BERKY, the EMTA assay platform detects the active form of the G α subunits rather than the surrogate measurement of G $\alpha\beta\gamma$ dissociation (*Galés et al., 2005; Masuho et al., 2015; Maziarz et al., 2020; Mende et al., 2018*), which can also detect non-productive binding as recently described for the V₂ engagement of G₁₂ (*Okashah et al., 2020*).

A potential caveat of EMTA is the use of common downstream effectors for all members of a given G protein family. Indeed, one cannot exclude that distinct members of a given family may display different relative affinities for their common effector. However, such differences are compensated by our data normalization that establishes the maximal response observed for a given subtype as the reference for this pathway (*Figure 3A*), as long as the number of the diversity of receptors included in the analysis is sufficient.

A second potential caveat of EMTA is that, when using heterologously expressed GPCRs and G proteins, some of the responses could result from favorable stoichiometries that may not exist under physiological conditions. It follows that such profiling represents the coupling possibilities of a given GPCR and not necessarily the coupling that will be observed in all cell types. Any couplings observed in such high-throughput studies requires further validation to conclude on their physiological relevance in cells or tissues of interest, and to form hypothesis for futures studies. Because we elected to use unmodified receptors (i.e.: not bearing any tags), the expression level of receptors could not be directly monitored. However, the double normalization method developed (see Materials and methods) allows quantitative comparison of coupling preferences across different receptors curtailing the influence of the assay response windows as well as receptor expression levels. Indeed, the double normalization allows ranking the coupling propensity of the receptors first as a function of the receptor which shows the strongest coupling to a specific G protein subtype, and then establishing the maximal response observed for a given G protein subtype as the reference for all G protein activated by a given receptor. In addition, as illustrated using the ET_A receptor as example, titrating receptor levels did not influence the pEC₅₀ for the activation of the different G protein coupled to this receptor (Figure 2-figure supplement 3B and Supplementary file 1C). Similarly, the pEC_{50} was not affected when titrating the amount of G protein subtype expressed (Figure 2—figure supplement 3A and Supplementary file 1B). As expected, only the amplitude of the response was affected.

It could be argued that overexpressing the G protein effectors (i.e.: p63-RhoGEF, Rap1GAP or PDZ-RhoGEF) used as sensors could influence the couplings observed. This potential caveat is mitigated by the fact that we used truncated part and/or modified versions of these effectors that limit the possibilities of interference with other components of the signaling machinery, and served essentially as a binding detector of the active forms of the G proteins (see Materials and methods). Supporting this notion, titrating the amount of the biosensor effector component did not affect the pEC₅₀ of G protein activation (*Figure 2—figure supplement 3C* and *Supplementary file 1D*).

Another limitation of the EMTA platform is the lack of a soluble effector protein selective for activated $G\alpha_s$ thus requiring tagging of the $G\alpha_s$ subunit (*Figure 1B*, bottom) and monitoring its dissociation from the plasma membrane. Yet, our data show that this translocation reflects G_s activation state, justifying its use in a G protein activation detection platform.

Finally, because EMTA is able to detect constitutive activity, high receptor expression levels may lead to an elevated basal signal level that may obscure an agonist-promoted response. Such an example can be appreciated for the A₁ receptor for which the agonist-promoted $G\alpha_{12}$ response did not reach the activation threshold criteria because of a very high constitutive activity level (*Figure 5A*). The impact of receptor expression on the constitutive activity and the narrowing on the agonist-promoted response is illustrated for $G\alpha_{q}$ activation by the 5-HT_{2C} (*Figure 5—figure supplement 1B*).

A limitation of any large-scale signaling study and drug discovery program is that ligands may elicit responses downstream of receptors other than the one under study. The development of a G_z/G_{15} quasi-universal biosensor enables efficient screening and detection of such polypharmacology and cross-talk. Using a combination of EMTA and appropriate pharmacological tools, we also proposed a systematic approach to distinguish off-target action of ligands from cross-talk. Interestingly, the cross-talk between the M_3 and CB receptors detected (*Figure 6*) may have physiological relevance since activation of muscarinic acetylcholine receptors has been shown to enhance the release of endocannabinoids in the hippocampus (*Kim et al., 2002*). The combined G_z/G_{15} biosensor should be particularly useful for early profiling of compound activity on safety panels and for the design of drugs displaying polypharmacology, an approach that is increasingly considered for the development of neuropsychiatric drugs (*Roth et al., 2004*).

The EMTA platform undoubtedly represents a novel tool-set that could be amenable for high throughput screening of small molecules and biologics across an array of signaling pathways, allowing for the discovery of functionally selective molecules or for GPCR deorphanization campaigns. The ability of the EMTA platform to quantitatively assess G protein coupling selectivity firmly expands the concept of functional selectivity and potential ligand bias beyond the dichotomic G protein vs. β arrestin view and provides plausible functional selectivity profiles that could be tested for their biological and pharmacological outcomes.

Materials and methods Cells

HEK293 clonal cell line (HEK293SL cells), hereafter referred as HEK293 cells, were a gift from S. Laporte (McGill University, Montreal, Quebec, Canada) and previously described (**Namkung et al., 2016**). HEK293 cells devoid of functional $G\alpha_s$ (ΔG_s), $G\alpha_{12}$ and $G\alpha_{13}$ ($\Delta G_{12/13}$), $G\alpha_q$, $G\alpha_{11}$, $G\alpha_{14}$ and $G\alpha_{15}$ ($\Delta G_{q/11}$) and, $G\alpha_i$, and $G\alpha_o$ ($\Delta G_{i/o}$) proteins were a gift from Dr. A. Inoue (Tohoku University, Sendai, Miyagi, Japan) and previously described (**Devost et al., 2017**; **Namkung et al., 2018**; **Schrage et al., 2015**; **Stallaert et al., 2017**). Cells were maintained in Dulbecco's Modified Eagle Medium (DMEM, Wisent, Saint-Jean-Baptiste, QC, Canada) supplemented with 10% fetal bovine serum (FBS, Wisent) and 1% antibiotics (100 U/mL penicillin and 100 µg/mL streptomycin (PS); Wisent). HL-60 cells were obtained from ATCC and maintained in RPMI 1640 medium containing L-Glutamine and 25 mM HEPES (Gibco) supplemented with 20% FBS (Wisent) and 1/100 volume PS (Wisent). Differentiation of HL-60 cells into neutrophil-like cells was induced by maintaining the cells in growth medium containing 1.3% DMSO (Bioshop) during 5 days. Cardiomyocytes derived from induced pluripotent stem cells (iPSCs; iCell Cardiomyocytes) were obtained from FUJIFILM Cellular Dynamics (Madison, WI, USA) and maintained in maintenance medium provided with the cells (special formulation by FujiFilm). Cells were grown at 37 °C in 5% CO₂ and 90% humidity and checked for mycoplasma contamination.

Plasmids and ebBRET biosensor constructs

Only human GPCRs and human G α subunits were used in this study. An open reading frame of each full-length GPCR was cloned into pcDNA3.1(+) expression plasmid. Except when otherwise specified, GPCRs sequences were devoid of epitope tags.

 $G\alpha_s$ -67-Rlucll (*Carr et al., 2014*), $G\alpha_{i1}$ -loop-Rlucll and GFP10-G γ_1 (*Armando et al., 2014*), $G\alpha_{i2}$ loop-Rlucll and βarrestin2-Rlucll (*Quoyer et al., 2013*), $G\alpha_{o8}$ -99-Rlucll (*Mende et al., 2018*), $G\alpha_q$ -118-Rlucll (*Breton et al., 2010*), $G\alpha_{12}$ -136-Rlucll and PKN-RBD-Rlucll (*Namkung et al., 2018*), $G\alpha_{13}$ -130-Rlucll (*Avet et al., 2020*), GFP10-G γ_2 (*Galés et al., 2006*), βarrestin1-Rlucll (*Zimmerman et al., 2012*), rGFP-CAAX (*Namkung et al., 2016*), EPAC (*Leduc et al., 2009*), MyrPB-Ezrin-Rlucll (*Leguay et al., 2021*), HA- β_2 AR (*Lavoie et al., 2002*), signal peptide-Flag-AT₁ (*Goupil et al., 2015*), and EAAC-1 (*Brabet et al., 1998*) were previously described. Full-length, untagged G α subunits, G β_1 and G γ_9 were purchased from cDNA Resource Center. GRK2 was generously provided by Dr. Antonio De Blasi (Istituto Neurologico Mediterraneo Neuromed, Pozzilli, Italy).

To selectively detect $G_{i/o}$ activation, a construct coding for aa 1–442 of Rap1 GTPase-activating protein (comprising a $G_{i/o}$ binding domain) fused to Rluc8, was sequence-optimized, synthetized and subcloned at TopGenetech (St-Laurent, QC, Canada). From this construct, a RlucII-tagged version of Rap1GAP (1-442) with a linker sequence (GSAGTGGRAIDIKLPAT) between Rap1GAP and RlucII was created by Gibson assembly in pCDNA3.1_Hygro (+) GFP10-RlucII, replacing GFP10. Three substitutions (i.e. S437A/S439A/S441A) were introduced into the Rap1GAP sequence by PCR-mediated mutagenesis. These putative (S437 and S439) and documented (S441) (*McAvoy et al., 2009*) protein kinase A phosphorylation sites were removed in order to eliminate any G_s-mediated Rap1GAP recruitment to the plasma-membrane.

To selectively detect $G_{q/11}$ activation, a construct encoding the G_q binding domain of the human p63 Rho guanine nucleotide exchange factor (p63RhoGEF; residues: 295–502) tagged with RlucII was done from IMAGE clones (OpenBiosystems; Burlington, ON, Canada) and subcloned by Gibson assembly in pCDNA3.1_Hygro (+) GFP10-RlucII, replacing GFP10. The G_q binding domain of p63RhoGEF and RlucII were separated by the peptidic linker ASGSAGTGGRAIDIKLPAT. N-term part containing palmitoylation sites maintaining p63 to plasma membrane and part of its DH domain involved in RhoA binding/activation (*Aittaleb et al., 2010; Aittaleb et al., 2011*) are absent of the sensor.

To selectively detect $G_{12/13}$ activation, a construct encoding the $G_{12/13}$ binding domain of the human PDZ-RhoGEF (residues: 281–483) tagged with RlucII was done by PCR amplification from IMAGE clones (OpenBiosystems) and subcloned by Gibson assembly in pCDNA3.1_Hygro (+) GFP10-RlucII, replacing GFP10. The peptidic linker GIRLREALKLPAT is present between RlucII and the $G_{12/13}$ binding domain of PDZ-RhoGEF. The sensor is lacking the PDZ domain of PDZ-RhoGEF involved in protein-protein interaction, as well as actin-binding domain and DH/PH domains involved in GEF activity and RhoA activation (*Aittaleb et al., 2010*).

The sequence of each EMTA biosensors is provided in the Supplementary file 5.

Transfection

For BRET experiments, HEK293 cells (1.2 mL at 3.5×10^5 cells per mL) were transfected with a fixed final amount of pre-mixed biosensor-encoding DNA (0.57 µg, adjusted with salmon sperm DNA; Invitrogen) and human receptor DNA. Transfections were performed using a polyethylenimine solution (PEI, 1 mg/mL; Polysciences, Warrington, PA, USA) diluted in NaCl (150 mM, pH 7.0; 3:1 PEI/ DNA ratio). Gelatin solution (1%; Sigma-Aldrich, Saint-Louis, Missouri) was used to stabilize DNA/PEI transfection mixes. Following addition of cells to the stabilized DNA/PEI transfection mix, cells were immediately seeded $(3.5 \times 10^4 \text{ cells/well})$ into 96-well white microplates (Greiner Bio-one; Monroe, NC, USA) and maintained in culture for the next 48 hr in DMEM containing 2% FBS and 1% PS. DMEM medium without L-glutamine (Wisent) was used for transfection of cells with mGluR to avoid receptor activation and desensitization. For Neutrophil-like differentiated HL-60 cells, cells were resuspended in electroporation medium (growth medium containing an extra 15 mM of HEPES pH 7.0) at 25 × 10^{6} cells/mL. Electroporation reactions were prepared by adding 50 μ L of DNA mastermix (20 μ g total of DNA adjusted with salmon sperm DNA, supplemented with 210 mM NaCl) to 200 µL of cell suspension and transferring into 0.4 cm gap electroporation cuvettes (Bio-Rad). The cells were electroporated at 350 µF/400 V using a Bio-Rad Gene Pulser II electroporation system, washed in electroporation medium, and seeded in 96-well plates at 0.8 × 10⁶ cells/well in 200 µL of growth medium. BRET assays were performed 6 hr post-electroporation. For iPSC Cardiomyocytes, cells were seeded in 96-well plates pretreated with fibronectin (10 μ g/ml 60 min; Sigma-Aldrich) at 3.5 \times 10⁴ cells /well. After 48 hr, attached iPSCs cells were transfected with the indicated biosensor components, using TransIT-LT1 reagent (Mirus; Madison, WI, USA), according to manufacturer recommendation. BRET assays were performed 48 hr after transfection.

For Ca²⁺ experiments, cells (3.5×10^4 cells/well) were co-transfected with the indicated receptor, with or without Ga₁₅ protein, using PEI and seeded in poly-ornithine-coated 96-well clear-bottom black microplates (Greiner Bio-one) and maintained in culture for the next 48 hr.

For BRET-based imagery, cells (4 × 10⁵ cells/dish) were seeded into 35 mm poly-d-lysine-coated glass-bottom culture dishes (Mattek Corporation; Ashland, MA, USA) in 2 ml of fresh medium and incubated at 37 °C in 5% CO₂, 3 day before imaging experiments. Twenty-four hours later, cells were transfected with EMTA ebBRET biosensors and the indicated receptor (i.e. p63-RhoGEF-RlucIl/rGFP-CAAX + $G\alpha_q$ and AT₁, Rap1GAP-RlucIl/rGFP-CAAX + $G\alpha_{i2}$ and D₂ or PDZ-RhoGEF-RlucIl/rGFP-CAAX + $G\alpha_{13}$ and TP α R) using X-tremeGENE 9 DNA transfection reagent (3:1 reagent/DNA ratio; Roche) diluted in OptiMEM (Gibco) and maintained in culture for the next 48 hr in DMEM containing 10% FBS and 1% PS.

Bioluminescence resonance energy transfer measurement

Enhanced bystander BRET (ebBRET) was used to monitor the activation of each $G\alpha$ protein, as well as β arrestin 1 and 2 recruitment to the plasma membrane. G α_s protein activation was measured between the plasma membrane marker rGFP-CAAX and human $G\alpha_s$ -Rlucll in the presence of human $G\beta_1$, $G\gamma_9$ and the tested receptor. $G\alpha_s$ downstream cAMP production was determined using the EPAC biosensor and GPBA receptor. $G\alpha_{i/a}$ protein family activation was followed using the selective-G_{i/a} effector Rap1GAP-RlucII and rGFP-CAAX along with the human $G\alpha_{i1}$, $G\alpha_{i2}$, $G\alpha_{oA}$, $G\alpha_{oB}$, or $G\alpha_{z}$ subunits and the tested receptor. $G\alpha_{a/11}$ protein family activation was determined using the selective- $G_{a/11}$ effector p63-RhoGEF-RlucII and rGFP-CAAX along with the human $G\alpha_{a}$, $G\alpha_{11}$, $G\alpha_{14}$, or $G\alpha_{15/16}$ subunits and the tested receptor. $G\alpha_{12/13}$ protein family activation was monitored using the selective- $G_{12/13}$ effector PDZ-RhoGEF-RlucII and rGFP-CAAX in the presence of either $G\alpha_{12}$ or $G\alpha_{13}$ and the tested receptor. The expression level of the G α subunits was monitored by western blot in HEK293 cells that endogenously expressed $G\alpha_{i1}$, $G\alpha_{i2}$, $G\alpha_{12}$, $G\alpha_{13}$, $G\alpha_{\alpha}$, $G\alpha_{11}$, $G\alpha_{14}$, and $G\alpha_s$ but not $G\alpha_{\alphaA}$, $G\alpha_{\alphaB}$, $G\alpha_{\alpha}$, $G\alpha_{$ and $G\alpha_{15}$ (Figure 2—figure supplement 6). $G\alpha_{12/13}$ -downstream activation of the Rho pathway was measured using PKN-RBD-RlucII or Ezrin-RlucII and rGFP-CAAX with the indicated receptor. Barrestin recruitment to the plasma membrane was determined using DNA mix containing rGFP-CAAX and Barrestin1-Rlucll with GRK2 or Barrestin2-Rlucll alone or with GRK2 and the tested receptor. Glutamate transporters EAAC-1 and EAAT-1 were systematically co-transfected with the mGluR to prevent receptor activation and desensitization by glutamate secreted in the medium by the cells (Brabet et al., 1998). All ligands were also tested for potential activation of endogenous receptors by transfecting the biosensors without receptor DNA. The G_{y}/G_{15} biosensor consists of a combination of the

following plasmids: rGFP-CAAX, Rap1GAP-RlucII, G α_2 , p63-RhoGEF-RlucII and G α_{15} . For G protein activation detection using the BRET-based G $\alpha\beta\gamma$ dissociation sensors, cells were co-transfected with untagged G β_1 and G α_q -118-RlucII, G α_{12} -136-RlucII or G α_{13} -130-RlucII with GFP10-G γ_1 , or G α_{11} -loop-RlucII, G α_{12} -loop-RlucII or G α_{08} -99-RlucII with GFP10-G γ_2 , along with the indicated receptor.

The day of the BRET experiment, cells were incubated in HBSS for 1 hr at room temperature (RT). Cells were then co-treated with increasing concentrations of ligand (see Appendix 1—key resources table and Supplementary file 2 for details) and the luciferase substrate coelenterazine prolume purple (1 µM, NanoLight Technologies; Pinetop, AZ, USA) for 10 min at RT. Plates were read on a Synergy Neo microplate reader (BioTek Instruments, Inc; Winooski, VT, USA) equipped with 410 ± 80 nm donor and 515 ± 30 nm acceptor filters or with a Spark microplate reader (Tecan; Männedorf, Switzerland) using the BRET² manufacturer settings. The BRET signal (BRET²) was determined by calculating the ratio of the light intensity emitted by the acceptor over the light intensity emitted by the donor. To validate the specificity of the biosensor responses, cells were pretreated in the absence or presence of either the $G\alpha_q$ inhibitor UBO-QIC (100 nM, 30 min; Institute for Pharmaceutical Biology of the University of Bonn, Germany), the $G\alpha_{i\alpha}$ inhibitor PTX (100 ng/mL, 18 hr; List Biological Laboratories, Campbell, California, USA) or the $G\alpha_s$ activator CTX (0–200 ng/mL, 4 hr; Sigma-Aldrich) before stimulation with agonist. For inverse agonist activity detection of A1 or CB1 receptors, cells were stimulated during 10 min with increasing concentrations of DPCPX or rimonabant, respectively. For ligand-cross receptor activation experiments, cells were pretreated for 10 min with increasing concentrations of antagonists or inverse agonist (eticlopride for D₂, WB4101 for $\alpha_{2A}AR$, atropine for muscarinic receptors and AM-630 for CB₁) before a 10 min stimulation with an EC₈₀ concentration of the indicated agonist. BRET was measured as described above. For the safety target panel ligand screen using the combined G_z/G_{15} sensor, basal ebBRET level was first measured 10 min following the addition of coelenterazine prolume purple (1 µM) and ebBRET level was measured again following a 10 min stimulation with a single dose of the indicated ligand (1 μ M for endothelin-1 and 10 μ M for all other ligands). Technical replicates for each receptor were included on the same 96-well plate. For kinetics experiment of ET_A activation, basal BRET was measured during 150 s before cells stimulation with either vehicle (DMSO) or 1 µM of endothelin-1 (at time 0 sec) and BRET signal was recorded each 30 s during 3570 s. For the validation of G12/13-mediated signal by new identified G12/13-coupled receptor using PKN- or Ezrin-based BRET biosensors, cells were pretreated or not with the $G\alpha_a$ inhibitor YM-254890 (1 µM, 30 min; Wako Pure Chemical Industries (Fujifilm), Osaka, Japan) before agonist stimulation for 10 min. For G protein activation detection using the BRET-based $G\alpha\beta\gamma$ dissociation sensors, and for titration experiments of either $G\alpha$ proteins subunit with GEMTA sensors, GPCRs with GEMTA sensors or Effector-Rlucll (p63-RhoGEF-Rlucll for $G\alpha_{\alpha/11}$, Rap1GAP-Rlucll for $G\alpha_{i/o}$ or PDZ-RhoGEF-Rlucll for $G\alpha_{12/13}$) from GEMTA sensors, cells were stimulated with increasing concentrations of the indicated agonist in the presence of prolume purple for 10 min before BRET measurement. For BRET in iPSC cardiomyocytes and HL-60 cells, cells were incubated in Tyrode Hepes buffer (137 mM NaCl, 0.9 mM KCl, 1 mM MgCl₂, 11.9 mM NaHCO₃, 3.6 mM NaH₂PO₄, 25 mM HEPES, 5.5 mM D-Glucose and 1 mM CaCl₂, pH 7.4) 30 min at RT before being treated with increasing concentrations of agonist for 15 min, using prolume purple (2 µM) as luciferase substrate, and BRET measured.

BRET data analyses and coupling efficiency evaluation

All BRET ratios were standardized using the equation below and represented as universal BRET (uBRET) values: uBRET = ((BRET ratio - A)/(B-A)) * 10,000. Constants A and B correspond to the following values:

A = pre-established BRET ratio obtained from transfection of negative control (vector coding for RlucII alone).

B = pre-established BRET ratio obtained from transfection of positive control (vector coding for a GFP10-RlucII fusion protein).

For a given signaling pathway, *u*BRET values at each agonist concentration were normalized as the % of the response obtained in the absence of agonist (vehicle) and concentration-response curves were fitted in GraphPad Prism 8.3 software using a four-parameter logistic nonlinear regression model. Results are expressed as mean ± SEM of at least three independent experiments.

A ligand-promoted response was considered real when the E_{max} value was \geq to the mean + 2*SD of the response obtained in vehicle condition and that a pEC₅₀ value could be determined in the agonist

concentration range used to stimulate the receptor. Consequently, a score of 0 or 1 was assigned to each signaling pathway depending on an agonist's ability to activate the tested pathway (0 = no activation; 1 = activation). In the case were responses associated to endogenous receptor were detectable, we considered as 'distorted' and exclude all the responses observed in the presence of transfected receptor for which E_{max} was \leq to 2*mean of the E_{max} value obtained with endogenous receptors or pEC₅₀ was \geq to 2*mean of the pEC₅₀ value obtained with endogenous receptors. Consequently, a score of 0 was assigned for these distorted responses in radial graph representation (*Figure 3—figure supplement 1*) and concentration-response curves were placed on a gray background in signaling signature profile panels (*Supplementary file 3*). Whenever transfected receptors, responses were considered 'true' and were assigned with a score of 1 for radial graph representation (*Figure 3—figure supplement 1*) and concentration-response curves were placed on a yellow background in signaling signature profile panels to indicate a partial contribution of endogenous receptors (*Supplementary file 3*).

We used a double normalization of E_{max} and pEC_{50} values to compare the signaling efficiency obtained for the 100 GPCRs across all receptors and pathways. E_{max} and pEC_{50} values deduced from concentration-response curves were first normalized between 0 and 1 across receptors by ranking the receptors as a function of the receptor that most efficiently activate a given pathway and then using the activation value for the pathway (including G protein and β arrestin subtypes) that a given receptor most efficiently activate as a reference for the other pathways that can be activated by this receptor. This double normalization can be translated in the following formalized equation:

• STEP1: For each receptor and for each pathway:

$$\begin{bmatrix} E_{max} GPCR_x \\ E_{max} GPCR_{Ref} \end{bmatrix}_{Pathway A} = Pathway specific normalized score for GPCR_x on pathway A ([PSNS GPCR_t]_{Pathway A})$$

where: $GPCR_x$ is receptor being analyzed, $GPCR_{Ref}$ is the receptor giving greatest E_{max} on pathway A of all receptors studied (i.e. reference receptor for pathway A). A PSNS was determined for every receptor and every pathway coupled to that receptor.

• STEP2: For any given receptor:

 $\frac{[PSNS GPCR_x]_{Pathway A}}{[PSNS GPCR_x]_{Pathway}} = \text{Normalized pathway A coupling score for GPCR}_x$

where: [PSNS GPCR_x] _{Pathway A} is the pathway specific normalized score for GPCR_x on pathway A, and [PSNS GPCR_x] _{Ref pathway} is the pathway specific normalized score for the pathway giving the highest PSNS for GPCR_x (i.e., reference pathway for GPCR_x).

For the safety target panel ligand screen using the combined G_z/G_{15} sensor, the fold ligand-induced stimulation was calculated for each receptor by dividing the BRET ratio after ligand addition (measured at 10 min post stimulation) by the basal BRET ratio prior to receptor stimulation. Activation thresholds were defined as the mean + 2*SD of the ligand-stimulated response obtained with receptor-null cells expressing only the combined G_z/G_{15} sensor.

Ca²⁺ mobilization assay

The day of the experiment, cells were incubated with 100 μ L of a Ca²⁺-sensitive dye-loading buffer (FLIPR calcium five assay kit, Molecular Devices; Sunnyvale, CA, USA) containing 2.5 mM probenecid (Sigma-Aldrich) for 1 hr at 37 °C in a 5% CO₂ incubator. During a data run, cells in individual wells were exposed to an EC₈₀ concentration of agonist, and fluorescent signals were recorded every 1.5 s for 3 min using the FlexStation II microplate reader (Molecular Devices). For receptors that also activate other G_{q/11} family members, cells were pretreated with the G_{q/11} inhibitor YM-254890 (1 μ M, 30 min) before agonist stimulation. G α_{15} is resistant to inhibition by YM-254890, thus allowing to measure Ca²⁺ responses generated specifically by G α_{15} .

BRET-based imaging

BRET images were obtained as previously described (*Kobayashi et al., 2019*). Briefly, the day of imaging experiment, cells were carefully rinsed with HBSS, and images were acquired before and after agonists addition (100 nM for Angiotensin II and U46619, and 1 μ M for dopamine) diluted in HBSS in the presence of the luciferase substrate coelenterazine prolume purple (20 μ M).

Images were recorded using an inverted microscope (Nikon Eclipse Ti-U) equipped with x60 objective lens (Nikon CFI Apochromat TIRF) and EM-CCD camera (Nuvu HNu 512). Measurements were carried out in photon counting mode with EM gain 3000. Exposure time of each photon counting was 100ms. Successive 100 frames were acquired alternatively with 480 nm longpass filter (acceptor frames) or without filter (total luminescence frames), and integrated. Image integrations were repeated 10 times and 1000 frames (*Video 1*) or 5 times and 500 frames (*Videos 2 and 3*) of acceptor and total luminescence were used to generate each image.

BRET values were obtained by dividing acceptor counts by total luminescence counts pixelwise. BRET values from 0.0 to 0.8 (*Video 1*) or 0.0–0.5 (*Videos 2 and 3*) were allocated to 'jet' heatmap array using MATLAB 2019b. Brightness of each pixel was mapped from the signal level of total luminescence image. 0% and 99.9% signal strength were allocated to the lowest and highest brightness to exclude the influence of defective pixels with gamma correction factor of 2.0.

The movies were generated using ImageJ 1.52 a. Frame rate is 10 (*Video 1*) or 3 (*Videos 2 and 3*) frames/s, and frame interval is 20 or 100 s for *Videos 1 and 2–3*, respectively. The field of view of the movie is 137 μ m x 137 μ m.

Western blot analysis

Cells were transfected or not with the indicated biosensors mix as previously described and whole-cell extracts were prepared 48 hr later. Briefly, cells were washed with ice-cold PBS and lysed in a buffer containing 10 mM Tris buffer (pH 7.4), 100 mM NaCl, 1 mM EDTA, 1 mM EGTA, 0.1% SDS, 1% Triton X-100, 10% Glycerol supplemented with protease inhibitors cocktails (Thermo Fisher Scientific). Cell lysates were centrifuged at 13,000 × g for 30 min at 4 °C. Equal amounts of proteins were separated by SDS-PAGE and transferred onto polyvinylidene fluoride membrane. The membranes were blocked (1 hr incubation at RT in TBS, 0.1% Tween-20, 5% BSA) and successively probed with primary antibody and appropriate goat secondary antibodies coupled to horseradish peroxidase (see **Appendix 1-Key Resources Table**). Western blots were visualized using enhanced chemiluminescence and detection was performed using a ChemiDoc MP Imaging System (BioRad). Relative densitometry analysis on protein bands was performed using MultiGauge software (Fujifilm). Results were normalized against control bands.

Statistical analyses

Curve fitting and statistical analyses were performed using GraphPad Prism 8.3 software and methods are described in the legends of the figures. Significance was determined as p < 0.05.

Acknowledgements

We thank Shane C Wright for scientific discussion and Monique Lagacé for critical reading of the manuscript. The authors are grateful to the funding from Bristol-Myers Squibb that supported the detection of $G\alpha$ proteins by endogenous receptors in iPSC cardiomyocytes and promyelocytic HL-60 cells.

Additional information

Competing interests

Arturo Mancini, Claire Normand, Florence Gross, Sandra Morissette: was employee of Domain Therapeutics North America during this research. Billy Breton: was employee of Domain Therapeutics North America during this research. Has filed patent application (US20200256869A1) related to the biosensors used in this work and the technology has been licensed to Domain Therapeutics. Christian Le Gouill, Hiroyuki Kobayashi, Mireille Hogue, Viktoriya Lukasheva: has filed patent application (US20200256869A1) related to the biosensors used in this work and the technology has been licensed to Domain Therapeutics. Eric B Fauman, Jean-Philippe Fortin: is employee and shares holders of Pfizer. Stephan Schann, Xavier Leroy: is employee and is part of the management of Domain Therapeutics. Michel Bouvier: is the president of Domain Therapeutics scientific advisory board. Has filed patent application (US20200256869A1) related to the biosensors used in this work and the technology has been licensed to Domain Therapeutics. The other authors declare that no competing interests exist.

Funding

Funder	Grant reference number	Author
Canada Research Chairs		Michel Bouvier
Canadian Institutes of Health Research	FDN-148431	Michel Bouvier
Lundbeckfonden	R218-2016-1266	David E Gloriam
Lundbeckfonden	R313-2019-526	David E Gloriam
Novo Nordisk Fonden	NNF18OC0031226	David E Gloriam
Lundbeckfonden	R278-2018-180	Alexander S Hauser
Bristol-Myers Squibb		Michel Bouvier

The funders had no role in study design, data collection and interpretation, or the decision to submit the work for publication.

Author contributions

Charlotte Avet, Conceptualization, Formal analysis, Investigation, Methodology, Writing – original draft, Writing – review and editing; Arturo Mancini, Conceptualization, Formal analysis, Investigation, Methodology, Resources, Writing – original draft, Writing – review and editing; Billy Breton, Conceptualization, Investigation, Methodology; Christian Le Gouill, Conceptualization, Methodology; Alexander S Hauser, Formal analysis, Writing – review and editing; Claire Normand, Hiroyuki Kobayashi, Florence Gross, Mireille Hogue, Viktoriya Lukasheva, Stéphane St-Onge, Marilyn Carrier, Sandra Morissette, Investigation; Madeleine Héroux, Investigation, Supervision; Eric B Fauman, Jean-Philippe Fortin, Resources; Stephan Schann, Funding acquisition, Resources; Xavier Leroy, Conceptualization, Resources, Supervision; David E Gloriam, Formal analysis, Funding acquisition, Supervision, Formal analysis, Funding acquisition, Methodology, Resources, Supervision, Methodology, Resources, Supervision, Methodology, Resources, Supervision, Writing – review and editing; Michel Bouvier, Conceptualization, Formal analysis, Funding acquisition, Methodology, Resources, Supervision, Writing – review and editing; Michel Bouvier, Conceptualization, Formal analysis, Funding acquisition, Methodology, Resources, Supervision, Writing – original draft, Writing – review and editing; Michel Bouvier, Conceptualization, Formal analysis, Funding acquisition, Methodology, Resources, Supervision, Writing – original draft, Writing – review and editing

Author ORCIDs

Alexander S Hauser (b http://orcid.org/0000-0003-1098-6419 Hiroyuki Kobayashi (b http://orcid.org/0000-0003-4965-0883 David E Gloriam (b http://orcid.org/0000-0002-4299-7561 Michel Bouvier (b http://orcid.org/0000-0003-1128-0100

Decision letter and Author response

Decision letter https://doi.org/10.7554/eLife.74101.sa1 Author response https://doi.org/10.7554/eLife.74101.sa2

Additional files

Supplementary files

Supplementary file 1. Impact of G proteins, receptor or effector-Rlucll titration on absolute pEC50 values. (A) Absolute pEC₅₀ values of responses elicited in WT vs. Knockout Gα protein background cells. pEC₅₀ values deduced from dose-response curves for various receptor tested in parental (WT) HEK293 cells or devoid of G_s, G_{12/13}, G_{q/11} or G_{i/o} proteins are related to *Figure 2—figure supplement* 1. (B) Absolute pEC₅₀ values deduced from dose-responses elicited in cells transfected with different amounts of Gα proteins. pEC₅₀ values deduced from dose-response curves obtained following Gα subunit titration in HEK293 cells transfected with GEMTA sensors and related to *Figure 2—figure supplement* 3A. (C) Absolute pEC₅₀ values of responses elicited in cells transfected with different amounts of receptors. pEC₅₀ values deduced from dose-response curves obtained following ET_A titration in HEK293 cells transfected with GEMTA sensors and related to *Figure 2—figure supplement* 3B. (D) Absolute pEC₅₀ values of responses elicited in cells transfected with different amounts of Effector-Rlucll. pEC₅₀ values deduced from dose-response curves obtained following

Effector-RlucII titration in HEK293 cells transfected with GEMTA sensors and related to **Figure 2 figure supplement 3C**.

• Supplementary file 2. List of tested receptors and ligands, along with the raw E_{max} , absolute pEC₅₀ and their corresponding double normalized (dnor) values. The E_{max} (in % of vehicle response) and absolute pEC₅₀ values deduced from concentration-response curves for the 100 GPCRs tested as well as the double normalized E_{max} and pEC₅₀ values calculated are related to **Supplementary file 3** and **Figure 3**, respectively.

• Supplementary file 3. Signaling profiles of 100 therapeutically-relevant human GPCRs using the EMTA ebBRET platform. Concentration-response curves in HEK293 cells expressing the indicated biosensor after stimulation of heterologously expressed receptor with the indicated ligand. Data are the mean ± SEM from at least 3 independent experiments and expressed in % of the response obtained for cells treated with vehicle. For ligands that elicited endogenous receptor-mediated responses (curves with light gray and yellow background for responses similar to and better responses than those obtained with the endogenous receptors, respectively), curves from cells expressing endogenous or heterologously expressed receptors are shown in *Figure 3—figure supplement 2*.

• Supplementary file 4. Comparison of G protein couplings identified with EMTA platform and other datasets. Comparison of G protein couplings identified with EMTA platform and TGF- α shedding assay in *Inoue et al., 2019* (A) or reported in GtP database (B).

- Supplementary file 5. EMTA biosensor amino acid sequences.
- Transparent reporting form

Data availability

All data generated or analysed during this study are included in the manuscript and supporting file; Source Data files have been provided for Figures 2, 4, 5, 6 and 7 and associated figure supplements. Supplementary file 1 contains the numerical data used to generate Figure 2-figure supplement 1 and Figure 2-figure supplement 3. Supplementary File 2 contains the numerical data used to generate figure 3. Further data are available in the companion paper co-published in eLife: https://doi.org/10.7554/eLife.74107.

References

- Aittaleb M, Boguth CA, Tesmer JJG. 2010. Structure and function of heterotrimeric G protein-regulated Rho guanine nucleotide exchange factors. *Molecular Pharmacology* **77**:111–125. DOI: https://doi.org/10.1124/mol. 109.061234, PMID: 19880753
- Aittaleb M, Nishimura A, Linder ME, Tesmer JJG. 2011. Plasma membrane association of p63 Rho guanine nucleotide exchange factor (p63RhoGEF) is mediated by palmitoylation and is required for basal activity in cells. The Journal of Biological Chemistry 286:34448–34456. DOI: https://doi.org/10.1074/jbc.M111.273342, PMID: 21832057
- Armando S, Quoyer J, Lukashova V, Maiga A, Percherancier Y, Heveker N, Pin J-P, Prézeau L, Bouvier M. 2014. The chemokine CXC4 and CC2 receptors form homo- and heterooligomers that can engage their signaling G-protein effectors and betaarrestin. FASEB Journal 28:4509–4523. DOI: https://doi.org/10.1096/fj.13-242446, PMID: 25053617
- Atwood BK, Lopez J, Wager-Miller J, Mackie K, Straiker A. 2011. Expression of G protein-coupled receptors and related proteins in HEK293, AtT20, BV2, and N18 cell lines as revealed by microarray analysis. BMC Genomics 12:14. DOI: https://doi.org/10.1186/1471-2164-12-14, PMID: 21214938
- Avet C, Sturino C, Grastilleur S, Gouill CL, Semache M, Gross F, Gendron L, Bennani Y, Mancini JA, Sayegh CE, Bouvier M. 2020. The PAR2 inhibitor I-287 selectively targets Galphaq and Galpha12/13 signaling and has anti-inflammatory effects. *Communications Biology* 3:719. DOI: https://doi.org/10.1038/s42003-020-01453-8
- Azzi M, Charest PG, Angers S, Rousseau G, Kohout T, Bouvier M, Piñeyro G. 2003. Beta-arrestin-mediated activation of MAPK by inverse agonists reveals distinct active conformations for G protein-coupled receptors. *PNAS* 100:11406–11411. DOI: https://doi.org/10.1073/pnas.1936664100, PMID: 13679574
- Bowes J, Brown AJ, Hamon J, Jarolimek W, Sridhar A, Waldron G, Whitebread S. 2012. Reducing safety-related drug attrition: the use of in vitro pharmacological profiling. *Nature Reviews. Drug Discovery* **11**:909–922. DOI: https://doi.org/10.1038/nrd3845, PMID: 23197038
- Brabet I, Parmentier ML, De Colle C, Bockaert J, Acher F, Pin JP. 1998. Comparative effect of L-CCG-I, DCG-IV and gamma-carboxy-L-glutamate on all cloned metabotropic glutamate receptor subtypes. *Neuropharmacology* 37:1043–1051. DOI: https://doi.org/10.1016/s0028-3908(98)00091-4, PMID: 9833633
- Breton B, Sauvageau É, Zhou J, Bonin H, Le Gouill C, Bouvier M. 2010. Multiplexing of multicolor bioluminescence resonance energy transfer. *Biophysical Journal* 99:4037–4046. DOI: https://doi.org/10.1016/j. bpj.2010.10.025, PMID: 21156147

- Bünemann M, Frank M, Lohse MJ. 2003. Gi protein activation in intact cells involves subunit rearrangement rather than dissociation. *PNAS* **100**:16077–16082. DOI: https://doi.org/10.1073/pnas.2536719100, PMID: 14673086
- Carr R, Du Y, Quoyer J, Panettieri RA, Janz JM, Bouvier M, Kobilka BK, Benovic JL. 2014. Development and characterization of pepducins as Gs-biased allosteric agonists. *The Journal of Biological Chemistry* 289:35668– 35684. DOI: https://doi.org/10.1074/jbc.M114.618819, PMID: 25395624
- Casey PJ, Fong HK, Simon MI, Gilman AG. 1990. Gz, a guanine nucleotide-binding protein with unique biochemical properties. The Journal of Biological Chemistry 265:2383–2390 PMID: 2105321.,
- Chandan NR, Abraham S, SenGupta S, Parent CA, Smrcka AV. 2021. Identification of G Protein α I Signaling Partners by Proximity Labeling Reveals a Network of Interactions That Includes PDZ-RhoGEF. *Pharmacology* and Toxicology. 1 :452545. DOI: https://doi.org/10.1101/2021.07.15.452545
- De Haan L, Hirst TR. 2004. Cholera toxin: a paradigm for multi-functional engagement of cellular mechanisms (Review)). Molecular Membrane Biology 21:77–92. DOI: https://doi.org/10.1080/09687680410001663267, PMID: 15204437
- Devost D, Sleno R, Petrin D, Zhang A, Shinjo Y, Okde R, Hebert TE. 2017. Conformational Profiling of the AT1 Angiotensin II Receptor Reflects Biased Agonism. G Protein Coupling, and Cellular Context. J Biol Chem 292:5443–5456. DOI: https://doi.org/10.1074/jbc.M116.763854
- Fukuhara S, Chikumi H, Gutkind JS. 2001. RGS-containing RhoGEFs: the missing link between transforming G proteins and Rho. Oncogene 20:1661–1668. DOI: https://doi.org/10.1038/sj.onc.1204182, PMID: 11313914
- Galandrin S, Oligny-Longpré G, Bouvier M. 2007. The evasive nature of drug efficacy: implications for drug discovery. Trends in Pharmacological Sciences 28:423–430. DOI: https://doi.org/10.1016/j.tips.2007.06.005, PMID: 17659355
- Galés C, Rebois RV, Hogue M, Trieu P, Breit A, Hébert TE, Bouvier M. 2005. Real-time monitoring of receptor and G-protein interactions in living cells. *Nature Methods* 2:177–184. DOI: https://doi.org/10.1038/nmeth743, PMID: 15782186
- Galés C, Van Durm JJJ, Schaak S, Pontier S, Percherancier Y, Audet M, Paris H, Bouvier M. 2006. Probing the activation-promoted structural rearrangements in preassembled receptor-G protein complexes. *Nature Structural & Molecular Biology* 13:778–786. DOI: https://doi.org/10.1038/nsmb1134, PMID: 16906158
- Goupil E, Fillion D, Clément S, Luo X, Devost D, Sleno R, Pétrin D, Saragovi HU, Thorin É, Laporte SA, Hébert TE. 2015. Angiotensin II type I and prostaglandin F2alpha receptors cooperatively modulate signaling in vascular smooth muscle cells. *The Journal of Biological Chemistry* **290**:3137–3148. DOI: https://doi.org/10. 1074/jbc.M114.631119, PMID: 25512374
- Hauser AS, Attwood MM, Rask-Andersen M, Schiöth HB, Gloriam DE. 2017. Trends in GPCR drug discovery: new agents, targets and indications. *Nature Reviews. Drug Discovery* 16:829–842. DOI: https://doi.org/10. 1038/nrd.2017.178, PMID: 29075003
- Hauser AS, Avet C, Normand C, Mancini A, Inoue A, Bouvier M, Gloriam DE. 2022. Common coupling map advances GPCR-G protein selectivity. *eLife* 11:e74107. DOI: https://doi.org/10.7554/eLife.74107, PMID: 35302494
- Hoffmann C, Gaietta G, Bünemann M, Adams SR, Oberdorff-Maass S, Behr B, Vilardaga J-P, Tsien RY, Ellisman MH, Lohse MJ. 2005. A FIAsH-based FRET approach to determine G protein-coupled receptor activation in living cells. Nature Methods 2:171–176. DOI: https://doi.org/10.1038/nmeth742, PMID: 15782185
- Inoue A, Raimondi F, Kadji FMN, Singh G, Kishi T, Uwamizu A, Ono Y, Shinjo Y, Ishida S, Arang N, Kawakami K, Gutkind JS, Aoki J, Russell RB. 2019. Illuminating G-Protein-Coupling Selectivity of GPCRs. *Cell* **177**:1933-1947. DOI: https://doi.org/10.1016/j.cell.2019.04.044, PMID: 31160049
- Jordan JD, Carey KD, Stork PJ, Iyengar R. 1999. Modulation of rap activity by direct interaction of Galpha(o) with Rap1 GTPase-activating protein. *The Journal of Biological Chemistry* **274**:21507–21510. DOI: https://doi.org/ 10.1074/jbc.274.31.21507, PMID: 10419452
- Kawamata Y, Fujii R, Hosoya M, Harada M, Yoshida H, Miwa M, Fukusumi S, Habata Y, Itoh T, Shintani Y, Hinuma S, Fujisawa Y, Fujino M. 2003. A G protein-coupled receptor responsive to bile acids. *The Journal of Biological Chemistry* 278:9435–9440. DOI: https://doi.org/10.1074/jbc.M209706200, PMID: 12524422
- Kenakin T. 2019. Biased Receptor Signaling in Drug Discovery. *Pharmacological Reviews* **71**:267–315. DOI: https://doi.org/10.1124/pr.118.016790, PMID: 30914442
- Kim J, Isokawa M, Ledent C, Alger BE. 2002. Activation of muscarinic acetylcholine receptors enhances the release of endogenous cannabinoids in the hippocampus. *The Journal of Neuroscience* **22**:10182–10191. DOI: https://doi.org/10.1523/JNEUROSCI.22-23-10182.2002, PMID: 12451119
- Kobayashi H, Picard L-P, Schönegge A-M, Bouvier M. 2019. Bioluminescence resonance energy transfer-based imaging of protein-protein interactions in living cells. *Nature Protocols* 14:1084–1107. DOI: https://doi.org/10. 1038/s41596-019-0129-7, PMID: 30911173
- Laschet C, Dupuis N, Hanson J. 2019. A dynamic and screening-compatible nanoluciferase-based complementation assay enables profiling of individual GPCR-G protein interactions. *The Journal of Biological Chemistry* 294:4079–4090. DOI: https://doi.org/10.1074/jbc.RA118.006231, PMID: 30593506
- Lavoie C, Mercier J-F, Salahpour A, Umapathy D, Breit A, Villeneuve L-R, Zhu W-Z, Xiao R-P, Lakatta EG, Bouvier M, Hébert TE. 2002. Beta 1/beta 2-adrenergic receptor heterodimerization regulates beta 2-adrenergic receptor internalization and ERK signaling efficacy. *The Journal of Biological Chemistry* 277:35402–35410. DOI: https://doi.org/10.1074/jbc.M204163200, PMID: 12140284

- Leduc M, Breton B, Galés C, Le Gouill C, Bouvier M, Chemtob S, Heveker N. 2009. Functional selectivity of natural and synthetic prostaglandin EP4 receptor ligands. *The Journal of Pharmacology and Experimental Therapeutics* 331:297–307. DOI: https://doi.org/10.1124/jpet.109.156398, PMID: 19584306
- Leguay K, Decelle B, He YY, Pagniez A, Hogue M, Kobayashi H, Le Gouill C, Bouvier M, Carréno S. 2021. Development of conformational BRET biosensors that monitor ezrin, radixin and moesin activation in real time. *Journal of Cell Science* **134**:jcs255307. DOI: https://doi.org/10.1242/jcs.255307, PMID: 33712451
- Lu M, Wang B, Zhang C, Zhuang X, Yuan M, Wang H, Li W, Su R, Li J. 2014. PQ-69, a novel and selective adenosine A1 receptor antagonist with inverse agonist activity. *Purinergic Signalling* **10**:619–629. DOI: https://doi.org/10.1007/s11302-014-9424-5, PMID: 25248972
- Lutz S, Shankaranarayanan A, Coco C, Ridilla M, Nance MR, Vettel C, Baltus D, Evelyn CR, Neubig RR, Wieland T, Tesmer JJG. 2007. Structure of Galphaq-p63RhoGEF-RhoA complex reveals a pathway for the activation of RhoA by GPCRs. Science (New York, N.Y.) 318:1923–1927. DOI: https://doi.org/10.1126/science.1147554, PMID: 18096806
- Mancini A, Frauli M, Breton B. 2015. Exploring the Technology Landscape of 7TMR Drug Signaling Profiling. Current Topics in Medicinal Chemistry 15:2528–2542. DOI: https://doi.org/10.2174/15680266156661507011 13344, PMID: 26126903
- Martin BR, Lambert NA. 2016. Activated G Protein Galphas Samples Multiple Endomembrane Compartments. The Journal of Biological Chemistry **291**:20295–20302. DOI: https://doi.org/10.1074/jbc.M116.729731, PMID: 27528603
- Masuho I, Ostrovskaya O, Kramer GM, Jones CD, Xie K, Martemyanov KA. 2015. Distinct profiles of functional discrimination among G proteins determine the actions of G protein-coupled receptors. *Science Signaling* 8:ra123. DOI: https://doi.org/10.1126/scisignal.aab4068, PMID: 26628681
- Maziarz M, Park JC, Leyme A, Marivin A, Garcia-Lopez A, Patel PP, Garcia-Marcos M. 2020. Revealing the Activity of Trimeric G-proteins in Live Cells with a Versatile Biosensor Design. *Cell* **182**:770-785.. DOI: https://doi.org/10.1016/j.cell.2020.06.020, PMID: 32634377
- McAvoy T, Zhou M, Greengard P, Nairn AC. 2009. Phosphorylation of Rap1GAP, a striatally enriched protein, by protein kinase A controls Rap1 activity and dendritic spine morphology. PNAS 106:3531–3536. DOI: https://doi.org/10.1073/pnas.0813263106, PMID: 19218462
- Mende F, Hundahl C, Plouffe B, Skov LJ, Sivertsen B, Madsen AN, Lückmann M, Diep TA, Offermanns S, Frimurer TM, Bouvier M, Holst B. 2018. Translating biased signaling in the ghrelin receptor system into differential in vivo functions. PNAS 115:E10255–E10264. DOI: https://doi.org/10.1073/pnas.1804003115, PMID: 30301804
- Meng J, Glick JL, Polakis P, Casey PJ. 1999. Functional interaction between Galpha(z) and Rap1GAP suggests a novel form of cellular cross-talk. The Journal of Biological Chemistry 274:36663–36669. DOI: https://doi.org/ 10.1074/jbc.274.51.36663, PMID: 10593970
- Namkung Y, Le Gouill C, Lukashova V, Kobayashi H, Hogue M, Khoury E, Song M, Bouvier M, Laporte SA. 2016. Monitoring G protein-coupled receptor and beta-arrestin trafficking in live cells using enhanced bystander BRET. Nature Communications 7:12178. DOI: https://doi.org/10.1038/ncomms12178, PMID: 27397672
- Namkung Y, LeGouill C, Kumar S, Cao Y, Teixeira LB, Lukasheva V, Giubilaro J, Simões SC, Longpré JM, Devost D, Hébert TE, Piñeyro G, Leduc R, Costa-Neto CM, Bouvier M, Laporte SA. 2018. Functional selectivity profiling of the angiotensin II type 1 receptor using pathway-wide BRET signaling sensors. *Science Signaling* 11:559. DOI: https://doi.org/10.1126/scisignal.aat1631
- Okashah N, Wright SC, Kawakami K, Mathiasen S, Zhou J, Lu S, Javitch JA, Inoue A, Bouvier M, Lambert NA. 2020. Agonist-induced formation of unproductive receptor-G12 complexes. *PNAS* **117**:21723–21730. DOI: https://doi.org/10.1073/pnas.2003787117, PMID: 32817560
- Oldham WM, Hamm HE. 2008. Heterotrimeric G protein activation by G-protein-coupled receptors. Nature Reviews. Molecular Cell Biology 9:60–71. DOI: https://doi.org/10.1038/nrm2299, PMID: 18043707
- Olsen RHJ, DiBerto JF, English JG, Glaudin AM, Krumm BE, Slocum ST, Che T, Gavin AC, McCorvy JD, Roth BL, Strachan RT. 2020. TRUPATH, an open-source biosensor platform for interrogating the GPCR transducerome. *Nature Chemical Biology* **16**:841–849. DOI: https://doi.org/10.1038/s41589-020-0535-8, PMID: 32367019
- Pándy-Szekeres G, Esguerra M, Hauser AS, Caroli J, Munk C, Pilger S, Keserű GM, Kooistra AJ, Gloriam DE. 2022. The G protein database, GproteinDb. *Nucleic Acids Research* 50:D518–D525. DOI: https://doi.org/10. 1093/nar/gkab852, PMID: 34570219
- **Quoyer J**, Janz JM, Luo J, Ren Y, Armando S, Lukashova V, Benovic JL, Carlson KE, Hunt SW, Bouvier M. 2013. Pepducin targeting the C-X-C chemokine receptor type 4 acts as a biased agonist favoring activation of the inhibitory G protein. *PNAS* **110**:E5088–E5097. DOI: https://doi.org/10.1073/pnas.1312515110, PMID: 24309376
- **Rojas RJ**, Yohe ME, Gershburg S, Kawano T, Kozasa T, Sondek J. 2007. Galphaq directly activates p63RhoGEF and Trio via a conserved extension of the Dbl homology-associated pleckstrin homology domain. *The Journal of Biological Chemistry* **282**:29201–29210. DOI: https://doi.org/10.1074/jbc.M703458200, PMID: 17606614
- Roth BL, Sheffler DJ, Kroeze WK. 2004. Magic shotguns versus magic bullets: selectively non-selective drugs for mood disorders and schizophrenia. Nature Reviews. Drug Discovery 3:353–359. DOI: https://doi.org/10.1038/ nrd1346, PMID: 15060530
- Sánchez-Soto M, Bonifazi A, Cai NS, Ellenberger MP, Newman AH, Ferré S, Yano H. 2016. Evidence for Noncanonical Neurotransmitter Activation: Norepinephrine as a Dopamine D2-Like Receptor Agonist. *Molecular Pharmacology* 89:457–466. DOI: https://doi.org/10.1124/mol.115.101808, PMID: 26843180

- Schrage R, Schmitz A-L, Gaffal E, Annala S, Kehraus S, Wenzel D, Büllesbach KM, Bald T, Inoue A, Shinjo Y, Galandrin S, Shridhar N, Hesse M, Grundmann M, Merten N, Charpentier TH, Martz M, Butcher AJ, Slodczyk T, Armando S, et al. 2015. The experimental power of FR900359 to study Gq-regulated biological processes. *Nature Communications* 6:10156. DOI: https://doi.org/10.1038/ncomms10156, PMID: 26658454
- Stallaert W, Westhuizen ET, Schonegge AM, Plouffe B, Hogue M, Lukashova V, Bouvier M. 2017. Purinergic Receptor Transactivation by the beta2-Adrenergic Receptor Increases Intracellular Ca(2+) in Nonexcitable Cells. *Molecular Pharmacology* 91:533–544. DOI: https://doi.org/10.1124/mol.116.106419
- Sunahara RK, Guan HC, O'Dowd BF, Seeman P, Laurier LG, Ng G, George SR, Torchia J, Van Tol HH, Niznik HB. 1991. Cloning of the gene for a human dopamine D5 receptor with higher affinity for dopamine than D1. *Nature* **350**:614–619. DOI: https://doi.org/10.1038/350614a0, PMID: 1826762
- Takasaki J, Saito T, Taniguchi M, Kawasaki T, Moritani Y, Hayashi K, Kobori M. 2004. A novel Galphaq/11selective inhibitor. *The Journal of Biological Chemistry* **279**:47438–47445. DOI: https://doi.org/10.1074/jbc. M408846200, PMID: 15339913
- Wedegaertner PB, Bourne HR, von Zastrow M. 1996. Activation-induced subcellular redistribution of Gs alpha. Molecular Biology of the Cell 7:1225–1233. DOI: https://doi.org/10.1091/mbc.7.8.1225, PMID: 8856666
- Wei H, Ahn S, Shenoy SK, Karnik SS, Hunyady L, Luttrell LM, Lefkowitz RJ. 2003. Independent beta-arrestin 2 and G protein-mediated pathways for angiotensin II activation of extracellular signal-regulated kinases 1 and 2. *PNAS* **100**:10782–10787. DOI: https://doi.org/10.1073/pnas.1834556100, PMID: 12949261
- Zimmerman B, Beautrait A, Aguila B, Charles R, Escher E, Claing A, Bouvier M, Laporte SA. 2012. Differential beta-arrestin-dependent conformational signaling and cellular responses revealed by angiotensin analogs. *Science Signaling* 5:ra33. DOI: https://doi.org/10.1126/scisignal.2002522, PMID: 22534132

Appendix 1.

Key Resources Table

Appendix 1—key resources table

Reagent type				
(species) or resource	Designation	Source or reference	Identifiers	Additional information
Cell line (Homo- sapiens)	HEK293	10.1038/ncomms12178 (Namkung et al., 2016)		HEK293 clonal cell line (HEK293SL cells)
Cell line (Homo- sapiens)	ΔG_s HEK293 cells	Dr. A. Inoue (Tohoku University, Sendai, Miyagi, Japan) 10.1124 / mol.116.106419 (Stallaert et al., 2017)		HEK293 cells devoid of functional Gα _s
Cell line (Homo- sapiens)	ΔG _{12/13} HEK293 cells	Dr. A. Inoue (Tohoku University, Sendai, Miyagi, Japan) 10.1074/jbc. M116.763854 (Devost et al., 2017)		HEK293 cells devoid of functional $G\alpha_{12}$ and $G\alpha_{13}$
Cell line (Homo- sapiens)	$\Delta G_{q^{r_{11}}}$ HEK293 cells	Dr. A. Inoue (Tohoku University, Sendai, Miyagi, Japan) 10.1038/ncomms10156 (<i>Schrage et al., 2015</i>)		HEK293 cells devoid of functional $G\alpha_{q_1}$ $G\alpha_{11}$, $G\alpha_{14}$ and $G\alpha_{15}$
Cell line (Homo- sapiens)	$\Delta G_{i/o}$ HEK293 cells	Dr. A. Inoue (Tohoku University, Sendai, Miyagi, Japan)		HEK293 cells devoid of functional Gα _i and Gα _o
Cell line (Homo- sapiens)	HL-60	ATCC	Cat. #: CCL-240	
Cell line (Homo- sapiens)	iCell Cardiomyocytes, 01434	FUJIFILM Cellular Dynamics	Cat. #: R1057	
Transfected construct (Homo sapiens)	Human Gα subunits-encoding plasmid library	Missouri S&T cDNA Resource Center (https://www.cdna.org/)	Cat. #: GNA1100000; GNA1200000; GNA0OA0000; GNA0OB0000; GNA0OB0000; GNA1200000; GNA1300001; GNA1400000; GNA1400000; GNA1400000; GNA1500000; GNA0SL0000	
Transfected construct (Homo sapiens)	Gβ1	Missouri S&T cDNA Resource Center (https://www.cdna.org/)	Cat. #: GNB0100000	
Transfected construct (Homo sapiens)	Gγ9	Missouri S&T cDNA Resource Center (https://www.cdna.org/)	Cat. #: GNG0900000	
Transfected construct (Homo sapiens)	Gαs-67-RlucII	10.1074/jbc.M114.618819(Carr et al., 2014)		
Transfected construct (Homo sapiens)	Gai1-loop-Rlucll	10.1096/fj.13–242446 (Armando et al., 2014)		
Transfected construct (Homo sapiens)	Gαi2-loop-Rlucll	10.1073/pnas.1312515110 (Quoyer et al., 2013)		
Transfected construct (Homo sapiens)	GαoB-99-RlucII	10.1073/pnas.1804003115(Mende et al., 2018)		

Reagent type (species) or resource	Designation	Source or reference	Identifiers	Additional information
Transfected construct (Homo	Cor 119 Bluell	10.1016 /j.bpj.2010.10.025 (Breton		
sapiens)	Gαq-118-Rlucll	et al., 2010)		
Transfected construct (Homo sapiens)	Gα12–136-RluclI	10.1126/scisignal.aat1631(Namkung et al., 2018)		
Transfected construct (Homo sapiens)	Ga13–130-Rlucll	10.1038 /s42003-020-01453-8 (Avet et al., 2020)		
Transfected construct (Homo sapiens)	GFP10-Gy1	10.1096/fj.13–242446 (Armando et al., 2014)		
Transfected construct (Homo				
sapiens)	GFP10-Gy2	10.1038/nsmb1134(Galés et al., 2006)		
Transfected construct (Homo sapiens)	EPAC	10.1124/jpet.109.156398(Leduc et al., 2009)		
Transfected construct (Homo sapiens)	rGFP-CAAX	10.1038/ncomms12178(Namkung et al., 2016)		
Transfected construct (Homo sapiens)	Rap1GAP-Rlucil	This paper		See Materials and Methods
Transfected construct (Homo				See Materials and Methods
sapiens)	p63-RhoGEF-RlucII	This paper		
Transfected construct (Homo sapiens)	PDZ-RhoGEF-RlucII	This paper		See Materials and Methods
Transfected construct (Homo sapiens)	PKN-RBD-RlucII	10.1126/scisignal.aat1631(Namkung et al., 2018)		
Transfected construct (Homo sapiens)	MyrPB-Ezrin-Rlucll	10.1242/jcs.255307(Leguay et al., 2021)		
Transfected construct (Homo		10.1126/scisignal.2002522(Zimmerman		
Transfected	βarrestin1-Rlucll	et al., 2012)		
(Homo sapiens)	βarrestin2-RlucII	10.1073/pnas.1312515110 (Quoyer et al., 2013)		
Transfected construct (Homo sapiens)	GRK2	This paper		See Materials and Methods
Transfected construct (Homo		· · · · · · · · · · · · · · · · · · ·		
sapiens)	RAMP3	Domain Therapeutics North America	N/A	
A 1.				

Reagent type (species) or resource	Designation	Source or reference	Identifiers	Additional information
Transfected construct (Homo		10.1016/s0028-3908(98)00091-4		
sapiens)	EAAC-1	(Brabet et al., 1998)		
Transfected construct (Homo sapiens)	EAAT-1	Domain Therapeutics North America	N/A	
			N/A	
Transfected construct (Homo sapiens)	signal peptide-Flag-AT ₁	10.1074/jbc.M114.631119(Goupil et al., 2015)		
Transfected construct (Homo				
sapiens)	FLAG-a2BAR	Domain Therapeutics North America	N/A	
Transfected construct (Homo sapiens)	HA-β ₂ AR	10.1074/jbc.M204163200(Lavoie et al., 2002)		
Antibody	Gαi1 (I-20) (Rabbit polyclonal)	Santa Cruz	Cat. #: sc-391 RRID: AB_2247692	WB (1:500)
Antibody	Gαi2 (T-19) (Rabbit polyclonal)	Santa Cruz	Cat. #: sc-7276 RRID:AB_2111472	WB (1:500)
Antibody	Gαo (K-20) (Rabbit polyclonal)	Santa Cruz	Cat. #: sc-387 RRID:AB_2111641	WB (1:500)
Antibody	Gaz (Rabbit monoclonal)	Abcam	Cat. #: ab154846	WB (1:1000)
Antibody	Gαs (K-20) (Rabbit polyclonal)	Santa Cruz	Cat. #: sc-823 RRID:AB_631538	WB (1:500)
Antibody	Gα12 (S-20) (Rabbit polyclonal)	Santa Cruz	Cat. #: sc-409 RRID:AB_2263416	WB (1:500)
Antibody	Gα13 (A-20) (Rabbit polyclonal)	Santa Cruz	Cat. #: sc-410 RRID:AB_2279044	WB (1:500)
Antibody	Gαq (E-17) (Rabbit polyclonal)	Santa Cruz	Cat. #: sc-393 RRID:AB_631536	WB (1:500)
Antibody	Gα11 (C-terminal) (Rabbit polyclonal)	Sigma-Aldrich	Cat. #: SAB2109181	WB (1:500)
Antibody	Gα14 (Rabbit polyclonal)	Sigma-Aldrich	Cat. #: SAB4300771	WB (1:500)
Antibody	Gα15 (Rabbit polyclonal)	ThermoFisher scientific (Pierce)	Cat. #: PA1-29022 RRID:AB_1958024	WB (1:5,000)
Antibody	βactin (Mouse monoclonal)	Sigma-Aldrich	Cat. #: A5441 RRID:AB_476744	WB (1:5,000)
Antibody	Anti-rabbit HRP-coupled (Donkey polyclonal)	GE Healthcare	Cat. #: NA934 RRID:AB_772206	WB (1:5,000)
Antibody	Anti-mouse HRP-coupled (Sheep polyclonal)	GE Healthcare	Cat. #: NA931 RRID:AB_772210	WB (1:10,000)
Commercial assay or kit	FLIPR Calcium 5 Assay Kit	Molecular Devices	Cat. #: R8185	
Chemical compound, drug	α-linolenic acid	Cayman Chemical	Cat. #: 21,910	
Chemical compound, drug	α-MSH	Sigma-Aldrich	Cat. #: M4135	
Chemical compound, drug	γ-MSH	Tocris	Cat. #: 4,272	
A			· · ·	

Reagent				
type (species) or resource	Designation	Source or reference	Identifiers	Additional information
Chemical compound,		. .	C + # 2.420	
drug	[Pyr1]-Apelin 13	Tocris	Cat. #: 2,420	
Chemical compound, drug	[Sar1, Ile4,8]-Angiotensin II	Peptides International	Cat. #: PAN-4476- V-1EA	
Chemical compound,				
drug	3-hydroxyoctanoic acid (3-HOA)	Sigma-Aldrich	Cat. #: H3898	
Chemical compound, drug	7α–25 dihydroxycholesterol	Sigma-Aldrich	Cat. #: SML0541	
Chemical				
compound, drug	Acetylcholine chloride	Tocris	Cat. #: 2,809	
Chemical compound, drug	ACT-389949	Provided by Bristol-Myers Squibb	N/A	
Chemical compound, drug	Adenosine	Sigma-Aldrich	Cat. #: A9251	
Chemical				
compound, drug	AM-630	Tocris	Cat. #: 1,120	
Chemical				
compound, drug	Amylin	Tocris	Cat. #: 3,418	
Chemical				
compound, drug	Angiotensin II (Ang II)	Sigma-Aldrich	Cat. #: A9525	
Chemical				
compound, drug	Arginine vasopressin (AVP)	Sigma-Aldrich	Cat. #: V9879	
Chemical compound,				
drug	Atropine	Sigma-Aldrich	Cat. #: A0132	
Chemical compound,				
drug	Bovine serum albumin	Sigma-Aldrich	Cat. #: A7030	
Chemical compound,				
drug	С5а	Complement Technology	Cat. #: A144(300)	
Chemical compound,			C	
drug Chemical	Calcitonin	Bachem	Cat. #: H-2250	
compound, drug	CCK Octapeptide, sulfated (CCK8)	Tocris	Cat. #: 1,166	
Chemical compound,				
drug	CCL20	R&D Systems	Cat. #: 360-MP/CF	
Chemical compound, drug	CCL3 (MIP-1a)	R&D Systems	Cat. #: 270-LD/CF	
Chemical				
compound, drug	Cholera Toxin (CTX) from Vibrio cholerae	Sigma-Aldrich	Cat. #: C8052	
Chemical compound, drug	Cmpd43	Provided by Bristol-Myers Squibb	N/A	
A	x 1 Continued on port page			

Appendi	x 1 Continued			
Reagent type (species) or resource	Designation	Source or reference	Identifiers	Additional information
Chemical compound, drug	Corticotropin-Releasing Factor (CRF)	Bachem	Cat. #: H-2435	
Chemical compound, drug	CXCL12	R&D Systems	Cat. #: 350-NS	
Chemical compound, drug	CXCL13	R&D Systems	Cat. #: 801 CX/CF	
Chemical compound, drug	CXCL8	R&D Systems	Cat. #: 208-IL/CF	
Chemical compound, drug	DAMGO	Tocris	Cat. #: 1,171	
Chemical compound, drug	Dopamine	Sigma-Aldrich	Cat. #: H8502	
Chemical compound, drug	DPCPX	Tocris	Cat. #: 0439	
Chemical compound, drug	Dynorphin A	Tocris	Cat. #: 3,195	
Chemical compound, drug	Endothelin-1	Tocris	Cat. #: 1,160	
Chemical compound, drug	Eticlopride	Tocris	Cat. #: 1,847	
Chemical compound, drug	Fingolimod	Provided by Bristol-Myers Squibb	N/A	
Chemical compound, drug	Gastric Inhibitory Peptide (GIP)	Bachem	Cat. #: H-5645	
Chemical compound, drug	Ghrelin	Tocris	Cat. #: 1,463	
Chemical compound, drug	Glucagon (Aittaleb et al., 2010; Aittaleb et al., 2011; Armando et al., 2014; Atwood et al., 2011; Avet et al., 2020; Azzi et al., 2003; Bowes et al., 2012; Brabet et al., 1998; Breton et al., 2010; Bünemann et al., 2003; Carr et al., 2014; Casey et al., 1990; Chandan et al., 2021; De Haan and Hirst, 2004; Devost et al., 2017; Fukuhara et al., 2001; Galandrin et al., 2007; Galés et al., 2005; Galés et al., 2006; Goupil et al., 2015; Hauser et al., 2017; Hauser et al., 2022; Hoffmann et al., 2005; Inoue et al., 2019; Jordan et al., 1999; Kawamata et al., 2003; Kenakin, 2019; Kim et al., 2002)	Bachem	Cat. #: H-6790	
Chemical compound, drug	Glucagon-like peptide-1 GLP-1 (Bowes et al., 2012; Brabet et al., 1998; Breton et al., 2010; Bünemann et al., 2003; Carr et al., 2014; Casey et al., 1990; Chandan et al., 2021; De Haan and Hirst, 2004; Devost et al., 2017; Fukuhara et al., 2001; Galandrin et al., 2007; Galés et al., 2005; Galés et al., 2006; Goupil et al., 2015; Hauser et al., 2017; Hauser et al., 2022; Hoffmann et al., 2005; Inoue et al., 2019; Jordan et al., 1999; Kawamata et al., 2003; Kenakin, 2019; Kim et al., 2002; Kobayashi et al., 2019; Laschet et al., 2019; Lavoie et al., 2002; Leduc et al., 2009; Leguay et al., 2021; Lu et al., 2014; Lutz et al., 2007)	Bachem	Cat. #: H-6795	

Reagent type (species) or resource	Designation	Source or reference	Identifiers	Additional information
Chemical compound, drug	Glucagon-like peptide-2 GLP-2 (Aittaleb et al., 2010; Aittaleb et al., 2011; Armando et al., 2014; Atwood et al., 2011; Avet et al., 2020; Azzi et al., 2003; Bowes et al., 2012; Brabet et al., 1998; Breton et al., 2010; Bimemann et al., 2003; Carr et al., 2014; Casey et al., 1990; Chandan et al., 2021; De Haan and Hirst, 2004; Devost et al., 2017; Fukuhara et al., 2001; Galandrin et al., 2017; Fukuhara et al., 2001; Gologi et al., 2015; Hauser et al., 2017; Hauser et al., 2022; Hoffmann et al., 2005; Inoue et al., 2019; Jordan et al., 1999; Kawamata et al., 2003; Kenakin, 2019; Kim et al., 2002; Kobayashi et al., 2019; Laschet et al., 2019; Lavoie et al., 2002; Leduc et al., 2009)	Bachem	Cat. #: H-7742	
Chemical compound, drug	Glutamate	Sigma-Aldrich	Cat. #: 49,621	
Chemical compound, drug Chemical	GnRH (LH-RH)	Peptides International	Cat. #: PLR-4013	
compound, drug	Histamine	Tocris	Cat. #: 3,545	
Chemical compound, drug	Kallidin	Anaspec	Cat. #: 22,853(AN)	
Chemical compound, drug	Leukotriene B4 (LTB4)	Cayman Chemical	Cat. #: 20,110	
Chemical compound, drug	Leukotriene D4 (LTD4)	Cayman Chemical	Cat. #: 20,310	
Chemical compound, drug	Litocholic acid	Sigma-Aldrich	Cat. #: L6250	
Chemical compound, drug	Melatonin	Bachem	Cat. #: Q-1300	
Chemical compound, drug	MDL 29,951	Cayman Chemical	Cat. #: 16,266	
Chemical compound, drug	Neuropeptide FF (NPFF)	Tocris	Cat. #: 3,137	
Chemical compound, drug	Neuropeptide Y (NPY)	Bachem	Cat. #: H-6375	
Chemical compound, drug	Nicotinic acid	Abcam	Cat. #: ab120145	
Chemical compound, drug	Nociceptin	Tocris	Cat. #: 910	
Chemical compound, drug	Noradrenaline	Tocris	Cat. #: 5,169	
Chemical compound, drug	Oleoyl-Lysophosphatidic acid (O-LPA)	Sigma-Aldrich	Cat. #: L7260	
Chemical compound, drug	Orexin-A	Bachem	Cat. #: H-4172	

Reagent type (species) or resource	Designation	Source or reference	Identifiers	Additional information
Chemical compound, drug	Oxytocin	Tocris	Cat. #: 1910	
Chemical	Parathyroid Hormone (Aittaleb et al., 2010; Aittaleb et al., 2011; Armando et al., 2014; Atwood et al., 2011; Avet et al., 2020; Azzi et al., 2003; Bowes et al., 2012; Brabet et al., 1998; Breton et al., 2010; Bünemann et al., 2003; Carr et al., 2014; Casey et al., 1990; Chandan et al., 2021; De Haan and Hirst, 2004; Devost et al., 2017; Fukuhara et al., 2001; Galandrin et al., 2007; Galés et al., 2005; Galés et al., 2006; Goupil et al., 2015; Hauser et al., 2017; Hauser et al., 2022; Hoffmann et al., 2005; Inoue et al., 2019; Jordan et al., 2005; Kawamata et al., 2003; Kenakin, 2019; Kim et al., 2005; Kobayashi et al., 2019; Laschet et al., 2019; Lavoie et al., 2002; Leduc et al., 2009; Leguay et al., 2021)	Sigma-Aldrich	Cat. #: P3796	
Chemical compound, drug	Pertussis toxin (PTX) from Bordetella pertussis	List Biological Laboratories	Cat. #: 179 A(LB)	
Chemical compound,				
Chemical compound, drug	pH (proton) (Hydrochloric acid) Probenecid	Sigma-Aldrich Sigma-Aldrich	Cat. #: 320,331 Cat. #: P8761	
Chemical compound, drug	Prolume Purple (methoxy e-Coelenterazine; Me-O-e- CTZ)	Nanolight	Cat. #: 369	
Chemical compound, drug	Propionate (sodium salt)	Sigma-Aldrich	Cat. #: P1880	
Chemical compound, drug	Prostaglandin D2 (PGD2)	Cayman Chemical	Cat. #: 12,010	
Chemical compound, drug	Prostaglandin E2 (PGE2)	Sigma-Aldrich	Cat. #: P0409	
Chemical compound, drug	RFamide-related peptide 3 (RFRP3)	Tocris	Cat. #: 4,683	
Chemical compound, drug	Rimonabant	Cayman Chemical	Cat. #: 9000484	
Chemical compound, drug	Saralasin	АрехВіо	Cat. #: B5063	
Chemical compound, drug	Serotonin	Cayman Chemical	Cat. #: 14,332	
Chemical compound, drug	SLIGKV-NH2 (PAR2 AP)	Tocris	Cat. #: 3,010	
Chemical compound, drug	SNC80	Sigma-Aldrich	Cat. #: S2812	
Chemical compound, drug	Somatostatin-14	Bachem	Cat. #: H-6276	
Chemical compound, drug	Sphingosine 1-phosphate	Cayman	Cat. #: 62,570	
A 1.				

properties concertsbackgrationSource or referenceidentifiantsAdditional informationampaired manuarityTrLRAND2 (PART AP)ToolsCat.#.1.444	Аррении	(Continued			
ampand trueTuckisTuckisCut. #: 1,644Jammad trugTKN027KN027KN027KN027Jennical ungIBC-OLIC (FR00339)Institute for Paramacutual Biology of trueN/AJennical ungIde-OLIC (FR00339)Institute for Paramacutual Biology of trueN/AJennical ungIde-OLIC (FR00339)Institute for Paramacutual Biology of trueK.A. #: 171.476Jennical ungIndecantic addSigma-AldrichCat. #: 171.476Jennical ungIncortin IISigma-AldrichCat. #: 171.476Jennical ungIncortin IISigma-AldrichCat. #: 171.476Jennical ungIncortin IISigma-AldrichCat. #: 171.476Jennical ungIncortin IISigma-AldrichCat. #: 171.476Jennical ungIncortin IITocrisCat. #: 171.476Jennical ungIncortin IITocrisCat. #: 171.476Jennical ungIndextro (Sigtim)Cat. #: 171.476Jennical ungIndextro (Sigtim)Cat. #: 171.476Jennical ungIndextro (Sigtim)Cat. #: 171.476Jennical ungIndextro (Sigtim)Cat. #: 171.476Jennical ungIndextro (Sigtim)Cat. #: 171.476Jennical ungIndextro (Sigtim)Cat. #: 172.40631Jennical ungIndextro (Sigtim)Cat. #: 171.476Jennical ungIndextro (Sigtim)Cat. #: 172.40631Jennical ungIndextro (Sigtim)Cat. #:	Reagent type (species) or resource	Designation	Source or reference	Identifiers	
ompoond transminuteFRV027Provided by Britet-Myers SqubbbN/ATammaid upBOC-GC (FR000359)Institute for Pharmacoutical Biology of the University of BonnN/ATammaid upUndecanoic acidSigma-AkrichCat # 171,476Semond ompoond upUncortin IIPhoenix PharmacouticalCat # 019-30Tammaid upUncortin IISigma-AkrichCat # 11066Tammaid upTarCisCat # 111Cat # 019-30Tammaid upVascactive Intestinal Peptide (MP)TorisCat # 111Tammaid 	Chemical compound, drug	TFLLR-NH2 (PAR1 AP)	Tocris	Cat. #: 1,464	
ompoor by Bo-Call C (FR00359)Institute for Pharmaceutical Biology the University of BonnNAPhennical impoord upon upon upon 	Chemical compound, drug	TRV027	Provided by Bristol-Myers Squibb	N/A	
ampound indecancic acidSigma AldrichCat. #: 171,776Demical homound ingenound ingenound ingenoundUndecancic acidPhoenix PharmaceuticalCat. #: 019-30Demical homound 	Chemical compound, drug	UBO-QIC (FR900359)		N/A	
ompound impound impoundUncortin IIPhoenix PharmaceuticalCat. #: 109–30hemical impound impoundCTPSigma-AldrichCat. #: 10106hemical impound rugVasoactive Intestinal Peptide (VIP)TocrisCat. #: 10106hemical 	Chemical compound, drug	Undecanoic acid	Sigma-Aldrich	Cat. #: 171,476	
ompound impound impound impound impound impound impound 	Chemical compound, drug	Urocortin II	Phoenix Pharmaceutical	Cat. #: 019–30	
ompound beneficial empound impoundVasocitive Intestinal Peptide (MP)ToorisCat. #: 1911beneficial ompound impound rugW4101ToorisCat. #: AdditionW1N55,212-2Enzo Life SciencesCat. #: BMLCR105beneficial 	Chemical compound, drug	UTP	Sigma-Aldrich	Cat. #: U1006	
ompord ingWB101TocrisCat #: 946ompord ompord(N55212-2)Enzo Life SciencesCat. #: BMLCR105behenical ompordvN254990(Wako Pure Chemical Industries) (Fujifin)Cat. #: 257-00631behenical ompordv264990Sigma-AldrichCat. #: 257-00631behenical ompordv2nc chloride (Zn2+)Sigma-AldrichCat. #: 259-00631oftware, 	Chemical compound, drug	Vasoactive Intestinal Peptide (VIP)	Tocris	Cat. #: 1911	
ompgond prigWINSS,212-2Enzo Life SciencesCat. #: BMLCR105ompgond prigWaSS,212-2Wako Pure Chemical Industries ("pillin")Cat. #: 257-00631ompgond prigond prigond prigond 	Chemical compound, drug	WB4101	Tocris	Cat. #: 946	
ompound IrugYM-254890Wake Pure Chemical Industries (Fujfilm)Cat. #: 257-00631ibminoud IrugCinc choride (Zn2+)Sigma-AldrichCat. #: 229,997oftware, IgorithPrisn, Version 8.3GraphPadCat. #: 229,997oftware, IgorithmMATLAB, Version R2019bMattWorksVersion 1000000000000000000000000000000000000	Chemical compound, drug	WIN55,212-2	Enzo Life Sciences	Cat. #: BMLCR105	
ompound lugZinc chloride (Zn2+)Sigma-AldrichCat. #: 229,997oftware, oftware, lgorithmMATLAB, Version 8.3GraphPadoftware, lgorithmMATLAB, Version R2019bMathWorksoftware, lgorithmImage J, Version 1.52 aNIH https://imagej.nih.gov/ij/oftware, lgorithmScipy, Version 1.52 aNIH https://www.scipy.orgoftware, 	Chemical compound, drug	YM-254890		Cat. #: 257–00631	
rugZinc chloride (Zn2+)Sigma-AldrichCat. #: 229,997oftware, IgorithmPrism, Version 8.3GraphPadoftware, IgorithmMATLAB, Version R2019bMathWorksoftware, IgorithmImageJ, Version 1.52 aNIH https://imagej.nih.gov/ij/oftware, IgorithmScipy, Version 1.4.1NIH https://imagej.nih.gov/ij/oftware, IgorithmScipy, Version 1.4.1Greiner Bio-oneCat. #: 655,083other96 W white plateGreiner Bio-oneCat. #: 655,090other96 W black plate, clear-bottomGreiner Bio-oneCat. #: 655,090other0ptiPlate-384, White Opaque 384-well MicroplatePerkin ElmerCat. #: 6007290other35 mm poly-d-lysine-coated glass-bottom culture dishesMattekCat. #: 405TSUSother10300e Digital DispenserTecanCat. #: 30097370otherSynergy NEO Luminescence microplate readerBioTek InstrumentsCat. #: 30097370otherFlexStation 2 Multi-mode, auto-pipetting microplatePerkin ElmerCat. #: 30097370otherInverted microscopeNikon CElipse Ti-UVerkerotherInverted microscopeNikon CElipse Ti-UVerker	Chemical				
Igorithm MATLAB, Version R2019b MathWorks Igorithm ImageJ, Version 1.52 a NIH https://imagej.nih.gov/ij/ Igorithm Scipy, Version 1.4.1 https://www.scipy.org oftware, Igorithm Scipy, Version 1.4.1 https://www.scipy.org Other 96 W white plate Greiner Bio-one Cat. #: 655,083 Other 96 W black plate, clear-bottom Greiner Bio-one Cat. #: 650,090 Other 96 W black plate, clear-bottom Greiner Bio-one Cat. #: 6007290 Other 96 W black plate, clear-bottom culture dishes Mattek Cat. #: 6007290 Other 35 mm poly-d-lysine-coated glass-bottom culture dishes Mattek Cat. #: 405TSUS Other J300e Digital Dispenser Tecan Cat. #: 30097370 Other Synergy NEO Luminescence microplate reader BioTek Instruments Cat. #: 30097370 Other FlexStation 2 Multi-mode, auto-pipetting microplate Molecular Devices Scient	drug	Zinc chloride (Zn2+)	Sigma-Aldrich	Cat. #: 229,997	
Igorithm ImageJ, Version 1.52 a NIH https://imagej.nih.gov/ij/ Igorithm Scipy, Version 1.4.1 https://www.scipy.org Igorithm Scipy, Version 1.4.1 Cat. #: 655,083 Other 96 W white plate Greiner Bio-one Cat. #: 655,083 Other 96 W black plate, clear-bottom Greiner Bio-one Cat. #: 655,090 Other 0ptiPlate-384, White Opaque 384-well Microplate Perkin Elmer Cat. #: 6007290 Other 35 mm poly-d-lysine-coated glass-bottom culture dishes Mattek P35GC-1.5-14 C Other 35 mm poly-d-lysine-coated glass-bottom culture dishes BioTek Instruments Cat. #: 405TSUS Other 78 + Dispensehead Cassettes Hp (Tecan) Cat. #: 30097370 Other 59,orgy, NEO Luminescence microplate reader BioTek Instruments Cat. #: 30097370 Other FlexStation 2 Multi-mode, auto-pipetting microplate Molecular Devices	Software, algorithm	Prism, Version 8.3	GraphPad		
Igorithm Scipy, Version 1.4.1 https://www.scipy.org Igorithm 96 W white plate Greiner Bio-one Cat. #: 655,083 Other 96 W black plate, clear-bottom Greiner Bio-one Cat. #: 655,090 Other 96 W black plate, clear-bottom Greiner Bio-one Cat. #: 6007290 Other OptiPlate-384, White Opaque 384-well Microplate Perkin Elmer Cat. #: 6007290 Other 35 mm poly-d-lysine-coated glass-bottom culture dishes Mattek P35GC-1.5–14 C Other Microplate washer BioTek Instruments Cat. #: 405TSUS Other D300e Digital Dispenser Tecan Tecan Other Synergy NEO Luminescence microplate reader BioTek Instruments Cat. #: 30097370 Other Synergy NEO Luminescence microplate reader BioTek Instruments Cat. #: 30097370 Other FlexStation 2 Multi-mode, auto-pipetting microplate molecular Devices Wolecular Devices Scient File Other Inverted microscope Nikon CFI Apochromat TIRF Ko objective lens Nikon CFI Apochromat TIRF	Software, algorithm	MATLAB, Version R2019b	MathWorks		
Igorithm 96 W white plate Greiner Bio-one Cat. #: 655,083 Other 96 W black plate, clear-bottom Greiner Bio-one Cat. #: 655,090 Other OptiPlate-384, White Opaque 384-well Microplate Perkin Elmer Cat. #: 6007290 Other 35 mm poly-d-lysine-coated glass-bottom culture dishes Mattek P35GC-1.5-14 C Other Microplate washer BioTek Instruments Cat. #: 405TSUS Other D300e Digital Dispenser Tecan Tecan Other Synergy NEO Luminescence microplate reader BioTek Instruments Cat. #: 30097370 Other FlexStation 2 Multi-mode, auto-pipetting microplate reader Molecular Devices View Cat. #: 405TSUS Other Inverted microscope Nikon CFI Apochromat TIRF View Cat. #: 405TSUS	Software, algorithm	ImageJ, Version 1.52 a	NIH https://imagej.nih.gov/ij/		
Dther96 W black plate, clear-bottomGreiner Bio-oneCat. #: 655,090DtherOptiPlate-384, White Opaque 384-well MicroplatePerkin ElmerCat. #: 6007290Dther35 mm poly-d-lysine-coated glass-bottom culture dishesMattekCat. #: P35GC-1.5-14 CDtherMicroplate washerBioTek InstrumentsCat. #: 405TSUSDtherD300e Digital DispenserTecanTecanDtherT8 + Dispensehead CassettesHp (Tecan)Cat. #: 30097370DtherSynergy NEO Luminescence microplate readerBioTek InstrumentsCat. #: 30097370DtherFlexStation 2 Multi-mode, auto-pipetting microplate readerMolecular DevicesStrumentsDtherInverted microscopeNikon CFI Apochromat TIRFKeo Objective lens	Software, algorithm	Scipy, Version 1.4.1	https://www.scipy.org		
OtherOptiPlate-384, White Opaque 384-well MicroplatePerkin ElmerCat. #: 6007290Other35 mm poly-d-lysine-coated glass-bottom culture dishesMattekCat. #: P35GC-1.5–14 COtherMicroplate washerBioTek InstrumentsCat. #: 405TSUSOtherD300e Digital DispenserTecanOtherT8 + Dispensehead CassettesHp (Tecan)Cat. #: 30097370OtherSynergy NEO Luminescence microplate readerBioTek InstrumentsCat. #: 30097370OtherFlexStation 2 Multi-mode, auto-pipetting microplate readerMolecular DevicesOtherInverted microscopeNikon Eclipse Ti-UOtherx60 objective lensNikon CFI Apochromat TIRF	Other	96 W white plate	Greiner Bio-one	Cat. #: 655,083	
Dther 35 mm poly-d-lysine-coated glass-bottom culture dishes Mattek Cat. #: P35GC-1.5-14 C Dther Microplate washer BioTek Instruments Cat. #: 405TSUS Dther D300e Digital Dispenser Tecan Dther T8 + Dispensehead Cassettes Hp (Tecan) Cat. #: 30097370 Dther Synergy NEO Luminescence microplate reader BioTek Instruments Cat. #: 30097370 Dther FlexStation 2 Multi-mode, auto-pipetting microplate reader Molecular Devices Cat. #: 30097370 Dther Inverted microscope Nikon Eclipse Ti-U Vikon Eclipse Ti-U Dther x60 objective lens Nikon CFI Apochromat TIRF Cat. #: 30097370	Other	96 W black plate, clear-bottom	Greiner Bio-one	Cat. #: 655,090	
35 mm poly-d-lysine-coated glass-bottom culture dishes Mattek P35GC-1.5-14 C Where Microplate washer BioTek Instruments Cat. #: 405TSUS Dther D300e Digital Dispenser Tecan Dther T8 + Dispensehead Cassettes Hp (Tecan) Cat. #: 30097370 Dther Synergy NEO Luminescence microplate reader BioTek Instruments Cat. #: 30097370 Dther FlexStation 2 Multi-mode, auto-pipetting microplate reader Molecular Devices Dther Inverted microscope Nikon CFI Apochromat TIRF	Other	OptiPlate-384, White Opaque 384-well Microplate	Perkin Elmer	Cat. #: 6007290	
Dther D300e Digital Dispenser Tecan Dther T8 + Dispensehead Cassettes Hp (Tecan) Cat. #: 30097370 Dther Synergy NEO Luminescence microplate reader BioTek Instruments Dther FlexStation 2 Multi-mode, auto-pipetting microplate reader Molecular Devices Dther Inverted microscope Nikon Eclipse Ti-U Dther x60 objective lens Nikon CFI Apochromat TIRF	Other	35 mm poly-d-lysine-coated glass-bottom culture dishes	Mattek		
Dther T8 + Dispensehead Cassettes Hp (Tecan) Cat. #: 30097370 Dther Synergy NEO Luminescence microplate reader BioTek Instruments Dther FlexStation 2 Multi-mode, auto-pipetting microplate reader Molecular Devices Dther Inverted microscope Nikon Eclipse Ti-U Dther x60 objective lens Nikon CFI Apochromat TIRF	Other	Microplate washer	BioTek Instruments	Cat. #: 405TSUS	
Other Synergy NEO Luminescence microplate reader BioTek Instruments Dther FlexStation 2 Multi-mode, auto-pipetting microplate reader Molecular Devices Dther Inverted microscope Nikon Eclipse Ti-U Dther x60 objective lens Nikon CFI Apochromat TIRF	Other	D300e Digital Dispenser	Tecan		
Dther FlexStation 2 Multi-mode, auto-pipetting microplate reader Molecular Devices Dther Inverted microscope Nikon Eclipse Ti-U Dther x60 objective lens Nikon CFI Apochromat TIRF	Other	T8 + Dispensehead Cassettes	Hp (Tecan)	Cat. #: 30097370	
reader Molecular Devices Dther Inverted microscope Nikon Eclipse Ti-U Dther x60 objective lens Nikon CFI Apochromat TIRF	Other	Synergy NEO Luminescence microplate reader	BioTek Instruments		
Other x60 objective lens Nikon CFI Apochromat TIRF	Other		Molecular Devices		
	Other	Inverted microscope	Nikon Eclipse Ti-U		
EMCCD camera Nuvu HNu 512	Other	x60 objective lens	Nikon CFI Apochromat TIRF		
	Other	EMCCD camera	Nuvu HNu 512		